Medio de cultivo Culture Medium

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For *Bacillus atrophaeus*, *Bacillus subtilis* and
Geobacillus stearothermophilus growth.

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
Parameters determined at time of manufacture according to ISO
11138-1: 2006, ISO 11138-2: 2006, ISO 11138-3: 2006, ISO 11138-4:
2006, 11138-5: 2006 and IRAM 37102-1,2 and 3: 1999.

ISO and USP Compliant
ATCC is registered trademark of American Type Culture Collection

Lic. Adrián J. Rovetto
Director Técnico
Quality Assurance Director

Industria Argentina - Made in Argentina

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English

Composition

Each screw capped bottle contains 500 ml of growth indicator
medium for *Bacillus atrophaeus* ATCC 9372, *Bacillus subtilis* DSM
5230 and *Geobacillus stearothermophilus* ATCC 7953.

Product description

MC1020-500 Bionova® culture medium is designed *Bacillus*
atrophaeus ATCC 9372, *Bacillus subtilis* DSM 5230 and *Geobacillus*
stearothermophilus ATCC 7953 growth.

If the sterilization process was not successful, indicator medium
will turn to yellow after incubation at 37±2 °C for 48 hours for
Bacillus atrophaeus ATCC 9372 spores subjected to Dry Heat and
EO processes; *Bacillus subtilis* DSM 5230 subjected to Steam at
Low Temperature; or incubation at 60±2 °C for 24 hours for
Geobacillus stearothermophilus ATCC 7953 spores subjected to
Steam/VH₂O₂ processes and 48 hours for Formaldehyde
processes, thus showing the presence of living *Bacillus*
atrophaeus ATCC 9372, *Bacillus subtilis* DSM 5230 or *Geobacillus*
stearothermophilus ATCC 7953 spores respectively.

If the sterilization process is successful indicator medium color
will remain the same. Final readout should be carried out after 48
or 24 hours of incubation at 37±2 °C or 60±2 °C, as appropriate.

Warning

MC1020-500 Bionova® culture medium should not be used for
other bacterial species growth different from *Bacillus atrophaeus*
ATCC 9372, *Bacillus subtilis* DSM 5230 or *Geobacillus*
stearothermophilus ATCC 7953.

Do not reuse.

Storage

Store it in a dark place at temperature between 4 °C and 20 °C,
30-80 % relative humidity.

Do not freeze.

Directions for use

1. Under sterile conditions transfer the Bionova® MC1020-500
culture medium to a sterile tube suitable for culturing bacteria.

2. Identify the tube containing Bionova® MC1020-500 culture
medium with a label indicating the sterilizer number (if there is
more than one), load number and processing date.

3. Open the envelope containing the spore carrier by cutting its
corner off, under sterile conditions (for example: laminar flow)
and transfer using a sterile forcep.

IMPORTANT: Use latex gloves and a mask when transferring the
MC1020-500 culture medium to the sterile tube as well as when
transferring the spores carrier from the envelope to the
MC1020-500 culture medium tube.

WARNING! Do not expose the carrier to any surface, including
the outside of the culture medium tube.

IMPORTANT: Everytime a processed carrier is incubated, a
non-sterilized spore carrier should be used as a positive control.
A positive control assures proper incubation conditions.

4. Incubate the spore carrier inside the MC1020-500 culture
medium at 37±2 °C or 60±2 °C as appropriate. A color change to
yellow of the growth indicator medium shows sterilization
process has failed. A final negative result is obtained if there is no
color change of the MC1020-500 culture medium after 48 or 24
hours appropriate, meaning sterilization process has been
successful. The color of the MC1020-500 indicator medium used
to incubate the non-sterilized spore carrier (positive control)
should turn to yellow for the result to be valid.

Register results and discard the tubes immediately, as it is
explained below.

WARNING! Do not reuse the sterilizer until the spores carrier
growth result is negative (MC1020-500 culture medium
containing a processed spores carrier remains the original color
after the 48 o 24 hours of incubation).

Disposal

Discard the culture medium tubes used for analyzing the spores
growth according to your country's healthcare and safety
regulations. MC1020-500 culture medium positive tubes can be
sterilized in gravity air displacement steam sterilizers at 121 °C
for at least 20 minutes or at 132 °C for at least 15 minutes, or at
134 °C for at least 10 minutes in vacuum-assisted steam
sterilizers.