
Indicadores Biológicos Biological Indicators

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
Biological Indicators For Formaldehyde sterilization

English

Composition

Each tube contains a population of *Geobacillus stearothermophilus* ATCC 7953 spores soaked on a carrier. It also has an growth indicator media of purple color contained in the glass ampoule.

Product description

Bionova® BT100 biological indicator is specifically designed for the monitoring of formaldehyde sterilization processes. If the sterilization process was not successful, the indicator media will change from purple to yellow after incubation at 60±2 °C, indicating the presence of live *Geobacillus stearothermophilus* spores. If the sterilization process was successful, the indicator media will remain purple after incubation. The final readout should be made after 48 hours of incubation at 60±2 °C.

Warning!

Do not use BT100 biological indicator for monitoring steam sterilization cycle, or ethylene oxide, dry heat, plasma or other sterilization processes.

Do not re-use the biological indicators.

Storage

Store it in a dark place at temperatures between 15 - 30 °C, 35 - 60 % relative humidity. Do not freeze.

Do not store these biological indicators near sterilizing agents or other chemical products.

Directions for use

1. Identify the Bionova® BT100 indicator by writing the sterilizer number (in case of having more than one sterilizer), load number, and processing date on the indicator label.
2. Pack the biological indicator along with materials to be sterilized in an appropriated package according to recommended sterilization practices. Place this package in those areas which a priori you consider most inaccessible for sterilizing agent.
3. Sterilize as usual.
4. After the sterilization process has finished, you can:
A. Remove the biological indicator from sterilization package for incubation prior to material aeration of the material.
B. Aerate the package containing the biological indicator with the rest of the load and then remove the biological indicator for incubation.
5. Check the chemical indicator on the label of the biological indicator. A color change to green confirms that the biological indicator has been exposed to Formaldehyde.

IMPORTANT: This color change does not indicate that the process was sufficient to achieve sterility. If the chemical indicator is unchanged, check the sterilization process.

6. Crush the glass ampule contained in the biological indicator and incubate at 60±2 °C.
- IMPORTANT:** Use a non-sterilized biological indicator as a positive control each time a processed biological indicator is incubated. The positive control ensure that correct incubation conditions were met.

7. Incubate the processed biological indicator and the indicator used as positive control for a maximum of 48 hours at 60±2 °C. Reading should be made at convenient intervals of 10 hours.

A color change from purple to yellow of the growth indicator media means a sterilization process failure has occurred. If after 48 hours there is no color change in the processed indicators, a final negative result is made (the sterilization process was acceptable). The positive control indicator should show a purple to yellow color change for the result to be valid.

NOTE: To extend the incubation time to more than 48 hours, we recommend using a humidified environment to prevent the complete evaporation of the culture medium contained in the indicator.

Record the positive ones and discard them immediately as it is shown below.

WARNING! Do not use the sterilizer until the biological indicator test results are negative (process indicator reads purple).

Disposal

Dispose of used biological indicators according to your country's healthcare and safety regulations. The positive biological indicator can be autoclaved at 121 °C for at least 20 minutes, or at 132 °C for 15 minutes in a gravity displacement steam sterilizer, or at 134 °C for 10 minutes in a vacuum assisted steam sterilizer.

Industria Argentina - Made in Argentina
Fabricado por: Terragene S.R.L. - Güemes 2879 - (2000) Rosario - Santa Fe - Argentina

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ATCC is registered trademark of American Type Culture Collection