## Multi-variable Indicators

### Results Reference Guide / Renk Tablosu

Unprocessed / Orjinal

Failures / Hatalar

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## Multi-variable Indicators CD20 For Steam sterilization processes

EN

### Type 4 indicators according to ISO 11140-1:2014

#### Product description

Chemdye® multi-variable indicators have been designed to react to Steam sterilization processes within loads, ensuring an adequate control of the effectiveness of sterilization processes (temperature, time and steam quality). Chemdye® CD20 multi-variable indicators have been manufactured with an innovative ink technology, highly sensitive to steam quality and air presence within the load.

The purple ink was developed to turn to green when the process reaches the stated values of the critical process variables. Multi-variable indicators allow to rapidly detect whether the critical parameters of the sterilization process for which they have been designed to respond have been reached.

Metals free

#### Characteristics

CD20: Chemical strips of 105 mm. x 18 mm.

Indicators 100% free of toxic heavy metals.

#### Indications for use

Use Chemdye® CD20 multi-variable indicators for monitoring steam sterilization processes: 121-134 °C.

#### Precautions

Do not expose Chemdye® multi-variable indicators to Ethylene Oxide, Dry Seat or any sterilization processes other than Steam.

#### Instructions for use


- Place a Chemdye® multi-variable indicator within each package, bag or tray that will be sterilized by steam.
- In rigid containers place the indicator in every corner, or at least two diagonally opposite corners, to ensure that the sterilization agent has reached the most inaccessible areas.
- After completing the process, analyze results. The indicator must turn to the reference color for considering that the indicated conditions were met.

#### Storage and shelf life

• Indicators should be stored in a dry place, protected from light at temperatures between 10-30 °C, 30-80 % Relative Humidity. Do not wet. Do not store close to sterilizing agents.

• Indicators have an expiration date of 5 years from the date of manufacture when stored under recommended conditions. Both expiration date and batch number of the product are specified on the package label and in each unit of the product. Do not use indicators after their expiration date.

• After use, the indicator will not change within 1 year when stored at previously indicated conditions.

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