Cleaning Indicators

Results reference guide

Unprocessed

Failures

Processed and correct

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Ultrasonic cleaning indicators

For ultrasonic cavitation performance tests

Control of cleaning and washing process of medical instruments and materials is highly important since the outcome affects subsequent disinfection and/or sterilization processes success. Ultrasonic washers are used to remove soil from joints, crevices, lumens, and other difficult to access locations. Cleaning takes place through a process called cavitation, tiny bubbles of vaporized liquid that rapidly collapse creating high pressure shock waves which dislodges soil from surfaces. These washers have special transducers that generate ultrasonic vibrations, necessary for creating the cavitation phenomenon. If any of these transducers is not functioning, the washer can have "cold spots", areas within the washer tank where there is no cavitation and, as consequence, ultrasonic cleaning becomes incomplete. Also, cavitation can be influenced by the presence of dissolved air on the tank water, temperature and hardness of tank water, use of cleaning solutions or detergents, among others, making a routine test of ultrasonic cleaning effectiveness indispensable.

Product description

Chemdyne® CDWU Indicator consists of a clear vial, with a reactive blue-colored solution and glass beads immersed in it. The formulation of the solution contained in the vial allows the monitoring of cavitation capacity of the ultrasonic washer. When cavitation functions properly, vibration of the glass pearls triggers a color change in the solution, from blue to yellow, through a range of green color intermediates. On the outside, the vial contains a transparent label with a printing of a cleaning indicator, allowing corroboration of CDWU indicator exposure to washing. This label, which also includes the lot number, was designed to detach easily from the tube and placed in the required documentation for the complete registration process. The CDWU indicator can be used for routine monitoring of ultrasonic cleaning processes and allows control of effectiveness of washing machines with different configurations wash.

Characteristics

CDWU indicator does not need a special holder to be used. Simply, place one vial on the ultrasonic machine zone you wish to monitor and proceed to perform the cleaning. CDWU indicator shows high stability against other similar products on the market, which ensures correct performance in the medium and long term, thereby guaranteeing a reliable result in the routine control of the ultrasonic cleaner.

CDWU indicator has high confidence and high reliability to reveal deficiencies in the ultrasonic cleaning process (cavitation).

The unique design of the CDWU indicator label not only allows for double-check of the washing process (through washing of the cleaning indicator ink impression), but also facilitates the elaboration of a complete documentation and registration. Washing of ink from the label indicates that the tube was exposed to the washing solution. Nevertheless, this change does not ensure the effectiveness of the process.

CDWU indicator can be used for monitoring ultrasonic cleaning processes in washing machines operating in a temperature range of 18 °C to 70 °C, with or without detergent or washing solution, and at 35 kHz or higher frequencies.

The materials used for manufacturing the CDWU indicator, such as vials and glass beads, can be washed and reused, making of this, a recyclable and economically profitable product.

Instructions for use

- 1-Take as many indicators as needed, according to the number of zones of the washing machine you want to monitor.
- 2-Shake the tube by inverting one time and use immediately.
- 3-Place the indicator(s) on the selected area.
- 4-Run the desired ultrasonic washing cycle.
- 5-Once the cycle is completed, remove the indicator(s) from the washing machine.
- 6-Analyze the results according to the result table and record the data obtained.

Warning: Degassing of the cleaning solution is extremely important to achieve satisfactory results. New or cooled solutions must be degassed before proceeding with the cleaning process. Degassing is performed after adding the detergent and is accomplished by turning on the ultrasound energy and/or increasing the temperature of the solution. Degassing is complete when small bubbles stop rising to the surface of the cleaning solution and you can see a pattern of ripples on the surface. The time required for degassing varies, depending on the capacity of the tank and temperature of solution, and can range from several minutes for a small tank to an hour or more for a large one (a tank without heating can take several hours to degas).

Results table

Not exposed Cavitation failures Correct cavitation

If unfavorable results are obtained, check:

- **Working frequency of the ultrasonic cleaner.** CDWU indicator can be used at frequencies higher than 35kHz. A color change slower than average will indicate a weak spot of cavitation in the cleaner tank. A negative result will point out a "cold spot" in that area of the tank.
- **Time and temperature of the ultrasonic washing.** CDWU indicator is designed to verify cavitation processes during ultrasonic washes with a minimum duration of 3 minutes (180 seconds), at temperatures between 18-70 °C.
- **Positioning and amount of wash load.** Excessive load (depending on size and ultrasonic cleaner manufacturer's recommendations), will reduce the capacity of the washing machine cavitation and, hence, of effectiveness of the washing process.
- **Presence of bubbles within the wash solution.** This is due to the presence of dissolved gases in the ultrasonic cleaner solution, which weakens the force of cavitation. We strongly recommend performing de-gassing of the cleaning solution before starting the ultrasonic cleaning.
- **Location of indicator.** It may indicate the presence of a "cold spot" or weak spot of cavitation in the ultrasonic washer.
- **General condition of ultrasonic washing machine.** Any visible damage in the washer, or unknown sound detected, could imply that the machine may not be working properly. After checking the actions above-mentioned, repeat the cleaning process. Regarding matters relative to washing machine operation, contact your technical service to check for possible causes of failure.

Observations

We strongly recommend the use of at least one CDWU indicator before each ultrasonic cleaning cycle.


The time it takes to the CDWU indicator to change its color from blue to yellow, is proportional to the power of cavitation of the ultrasonic washer; frequencies above 40kHz will cause a faster color change.

Ultrasonic washing machines should be checked at least once a week, and the results recorded. The control tests fall into two categories depending on the test in question, and the distribution of the indicators will be different according to the test to be performed (see figure). Both types of test are carried out without load. The periodic functional test is performed for the initial set-up of ultrasonic washing machine, and then quarterly or after repairs. The objective of this test is to verify the uniform distribution of cavitation in the washing machine and also helps identify "cold spots" in the tank. The routine test, however, is done weekly to verify that cavitation is taking place properly. This test will detect a possible deterioration in the performance of the ultrasonic cleaner and help to prevent ineffective washing.

Periodic functional test

Routine test

Precautions

Do not re-use the CDWU indicator.

Do not open the indicator before its use, since this will cause its deterioration.

Do not use the CDWU indicator for monitoring normal washing using non-ultrasonic cleaners.

Use gloves to wash and recycle CDWU indicator material, after being used in an ultrasonic washing process (vials and glass beads).

Storage

CDWU indicators should be stored at temperatures between 10-30 °C, specially protected from light and moisture, within its original packaging.