Liquid Peracetic Acid

Quality Certification

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Product description

Integron® IT400 Indicators are designed to react in High Level Disinfection processes that use Peracetic Acid (PAA) as disinfecting agent, assuring an adequate control of the disinfection process efficacy (temperature, time and PAA concentration). The reactive ink was developed for processes that work at temperatures of 50-60 °C and turns to fuchsia when Peracetic Acid concentration and time exposure are within the values shown in the certificate of analysis.

Characteristics

IT400 Chemical Indicators consist of an indicator strip of 10 mm x 45 mm, printed with a purple reactive ink specific for PAA, contained within an envelope of non absorbent material of 25 mm x 70 mm, which avoids direct contact of the reactive ink with the aqueous medium.

Metals free

100% free of all metals indicators

Indications for use

Use Integron® IT400 Chemical Indicators for monitoring Liquid Peracetic High Level Disinfection cycles.

Precautions

Do not expose this product to Steam, EO, Dry Heat, or any sterilization or disinfection processes other than Liquid PAA.

ATTENTION! Processes at temperatures below 50 °C might affect diffusion and sensing of the indicators, bringing on incorrect results.

Instructions of use

1. With the reusable silicon strip provided along with the product, place the envelope with the chemical test strip inside the equipment, together with the load to be processed with PAA. **WARNING!** The indicator must be completely submerged in the washing solution to work properly.

2. Process the load according to the established procedures.

3. After completing the process, analyze the results. The indicator must change to the reference color for considering that the desired disinfection conditions were met.

Storage and shelf life

• Keep this indicator in dry places, protected from light and at temperatures between 10-30 °C and 30-80 % relative humidity. Do not wet. Do not storage close to sterilizing or disinfection agents.

• Indicators have an expiration date of 5 years from the date of manufacture when stored in the recommended conditions. Both expiration date and batch of the product are specified on the box label and on each unit of product. Do not use indicators after their expiration date.

• After use, the indicator will not change in 1 year, if stored at the temperature and relative humidity previously indicated.

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