Indicador Multi-variable Multi-variable Indicator

Results Reference Guide

Unprocessed

Failures

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Multi-variable Indicator CD16

For Ethylene Oxide sterilization processes

Type 4 Indicator according to ISO 11140-1:2014

Product description

Chemdye® Multi-variable Indicators have been designed to react to Ethylene Oxide sterilization processes within loads, ensuring an adequate control of the effectiveness of sterilization processes (temperature, time, relative humidity and gas concentration). Purple/Brown ink was developed to turn to green when the process reaches the stated values of the critical process variables. Multi-variable Indicators allow to rapidly detect whether the critical parameters of the sterilization process for which they have been designed to respond have been reached.

Characteristics

CD16: Punched chemical strips of 200 x 15 mm, this allows use of the entire strip or half strip depending on the size of the package to be sterilized.

Metals free

Indicators 100% free of toxic heavy metals.

Indications for use

Use Chemdye® Multi-variable Indicators for monitoring Ethylene Oxide sterilization processes. Parameters specifications: 600 mg/l EO, RH 60%, 30 °C, 30 minutes // 600 mg/l EO, RH 60%, 54 °C, 20 minutes.

Precautions

Do not expose Chemdye® Multi-variable Indicators to Steam, Dry heat or any sterilization processes other than Ethylene Oxide.

Instructions for use

- Place a Chemdye® Multi-variable Indicator within each package, bag or tray that will be sterilized by Ethylene Oxide.
- In rigid containers place the indicator in every corner, or at least two diagonally opposite corners, to ensure that the sterilization agent has reached the most inaccessible areas.
- After completing the process analyze the results. The indicator must turn to the reference color for considering that the indicated conditions were reached.

Storage and shelf life

- Indicators should be stored in a dry place, protected from light at temperatures between 10-30 °C, 30-80 % relative humidity. Do not wet. Do not store near sterilizing agents.
- Indicators have an expiration date of 5 years from the date of manufacture when stored in recommended conditions. Both the expiration date and batch number of the product are specified on the package label and in each unit of the product. Do not use indicators after their expiration date.
- After use, the indicator will not change within 1 year when stored at previously indicated conditions.

Indicador Multi-variable Multi-variable Indicator

Results Reference Guide


Unprocessed

Failures

One or more sterilization parameters were not achieved. The materials cannot be considered sterile and must be processed again.

The main failure causes are:

- Inadequate cycle duration
- Improper cycle selection
- Overloading / incorrect loading of the sterilizer
- Improper humidity
- Inadequate sterilization agent concentration.

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