## Spores Suspension

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Survival time = (log<sub>10</sub> labeled population - 2) x labeled D-value

Tiempo de muerte / \_\_\_\_\_ seg.  
Kill time

Kill time = (log<sub>10</sub> labeled population + 4) x labeled D-value

Performance values for Steam and Formaldehyde were calculated using a filter paper inoculated with 0.1 ml of spore suspension as a BI. For VH<sub>2</sub>O<sub>2</sub>, a polypropylene tube inoculated with 0.1 ml of spore suspension was used as a BI.

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Parameters determined at time of manufacture according to ISO 11138-1: 2006, ISO 11138-3: 2006, ISO 11138-5: 2006 and IRAM 37102: 1999 (Parts 1 and 3) standards. The values shown are reproducible only under the same conditions under which they were determined.

ISO and USP Compliant

ATCC is a registered trademark of American Type Culture Collection

Lic. Adrian J. Rovetto  
Director Técnico  
Technical Director

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*Geobacillus stearothermophilus* ATCC 7953

## Composition

Each Bionova® BT20S Spore Suspension BT20S contains 10ml of bacterial spores of *Geobacillus stearothermophilus* ATCC 7953 in 40% of ethanol.

## Product description

Bionova® BT20S Spore Suspension has been designed for direct inoculation of products or for the preparation of customized Biological Indicators used for control of Steam, Formaldehyde, Hydrogen Peroxide and/or Peracetic Acid sterilization processes. The combined use of inoculated products, pre-sterilization microorganism count determinations (bioburden), and physical process measurements and controls should be employed to validate sterilization processes.

**WARNING:** Use BT20S biological indicator only for the control of Steam, Formaldehyde, Hydrogen Peroxide and/or Peracetic Acid sterilization processes. Do not use for other processes such as Ethylene Oxide, Dry Heat and/or Radiation sterilization.

## Storage

Store in a dark place and at temperatures between 2-8 °C, 35-60 % Relative Humidity. Do not freeze.

Do not store near sterilizing agents or other chemical products.

## Directions for use

- Perform inoculation operations in a clean area.
- Samples to be inoculated should be representative of product being sterilized.
- Select the spore suspension of appropriate population for the sterilization method use. For most purposes, inoculation of product with approximately one million (10<sup>6</sup>) spores provides a suitable challenge.
- NOTE: Spore suspensions are standardized on the basis of number of spores per 0.1 mL of suspension.
- Vortex vial before each use.
- Use a sterile pipette or syringe to accurately deliver the volume of suspension to be utilized.
- If syringe is used, disinfect stopper surface and pull syringe plunger halfway back. Insert needle through stopper, push plunger in, and then slowly withdraw plunger to fill syringe to desired volume.
- If a pipette is used, remove stopper and insert pipette. Dispense the desired volume.
- Inoculate the product appropriately. The area to be inoculated should be the most difficult for the sterilizing agent to penetrate.
- Identify the inoculated samples.
- Distribute "Inoculated Test Samples" throughout the sterilizer load.
- After sterilization cycle is completed, test the inoculated products as soon as possible using an appropriate culture medium, and incubate at 60 °C. We suggest Bionova® MC20 or MC1020 media to final reading at 24 hours for Steam/Hydrogen Peroxide sterilization and at 48 hours for Formaldehyde sterilization. For conventional culture medium like TSB, incubation for 7 days at 60 °C is recommended.

**IMPORTANT:** Each time a processed product is incubated, use an unprocessed product sample as a positive control. The positive control ensures that the incubation conditions were met. An uninoculated product can be used as a negative control. The positive vials have to be registered and then have to be processed as indicated below.

**WARNING:** Do not use the sterilizer until the biological indicator shows a negative result.

## Disposal

Discard Bionova® BT20S spore suspension and growth medium tubes according to your country's healthcare and safety regulations. Positive growth medium tubes can be sterilized in gravity air displacement steam sterilizers at 121 °C for at least 20 minutes or at 132 °C for at least 15 minutes, or at 134 °C for at least 10 minutes in vacuum-assisted steam sterilizers.