Spore ampoules

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For Steam sterilization

Product description

Bionova® BT22 Biological Indicator has been designed for monitoring Steam sterilization processes at temperatures between 121°C and 134°C. The device consists of a glass ampoule containing spores of the non-pathogenic microorganism *Geobacillus stearothermophilus* ATCC 7953. The culture medium has been specifically designed to change its color when microbial growth occurs.

The spore population contained in the biological indicator is killed after an efficient steam sterilization cycle. At lower temperatures or lower exposure times, some spores may survive and thus might be capable of growing. This growth can be detected by the indicator medium which turns from violet to yellow after 48 hours of incubation at 60°C. Conversely, if the sterilization process has been effective, there is no microbial growth and the ampoule medium maintains its original color.

Directions for use

1. Take the Bionova® BT22 Ampoules out of the refrigerator and leave them stabilize at room temperature.

2. Place the Bionova® BT22 Ampoules inside a container with the liquid being sterilized or inside a similar container with liquid. This container should be the same size and have the same volume as the one containing the liquid being sterilized.

3. Sterilize as usual.

4. Once the sterilization process has finished, open the sterilizer door, wait five minutes and remove the biological indicator from the container.

WARNING! Wear safety glasses and gloves when removing the BT22 Biological Indicator from the sterilized container.

5. Incubate the processed biological indicator, a negative control and a non-sterilized biological indicator as a positive control, at 60±2 °C. The positive control ensures that the incubation conditions were met.

The negative control, not being exposed to the sterilization process, can be used as a reference color to compare positive and negative results. Conversely, when used with the ampoules during the process, it ensures that any color change produced immediately after sterilization is the result of the process and not due to germination and/or spore growth.

Readout should be performed at convenient intervals of 12 hours.

It is recommended the use of Dual Bionova® IC10/20 or similar devices that assure an efficient temperature control.

WARNING! Traditional culture stoves have thermal profiles that differ from the above recommended incubators, thus their use might affect or delay the Bionova® BT22 Biological Indicators result.

Interpretation

If the sterilization process is adequate, *G. stearothermophilus* spores do not survive and after incubation the content of the ampoule remains violet. Conversely, if sterilization is not efficient the spores survive and the indicator medium turns to yellow indicating a failure in the process.

After 24 hours, a yellow-pink intermediate color might be observed depending on the integrity of the spores. The content of the vial can become turbid or the appearance of a precipitate can be observed due to the bacterial growth.

Final readout must be performed after 48 hours of incubation at 60±2 °C.

For results to be valid, the positive control indicator should turn to yellow and the negative control should not change its original color.

Record the positive results and discard them immediately as it is shown below.

WARNING! Do not re-use the sterilizer until biological indicator test results are negative (processed indicator culture medium remains its original color).

Storage

The ampoules should be stored away from light, refrigerated at temperatures between 4 - 8 °C and 30-80% Relative Humidity.

Storage at room temperature (up to 25 °C) is possible for a period no longer than 15 days.

Do not freeze.

Do not store close to other sterilization agents or chemical products.

Important: If the product is not stored at recommended temperatures the stability of the vial content might be affected.

WARNING!

Do not use Bionova® BT22 Biological Indicators for monitoring Dry Heat, Radiation, Ethylene Oxide, Chemical Vapors or any other sterilization process different from Steam.

The ampoule is fragile, protect it from shocks, falls and sudden temperature changes.

Check the color of the indicator before starting a sterilization process. It must be violet.

Do not use an indicator after the expiration date indicated in the package.

Each indicator is intended for a single use and must be discarded after use.

The indicator must be processed by a trained healthcare professional.

Disposal

Discard Bionova® BT22 Biological Indicators after use according to your country's healthcare and safety regulations. The positive biological indicator can be sterilized in gravity air displacement steam sterilizers at 121 °C for at least 20 minutes or at 132 °C for at least 15 minutes, or at 134 °C for at least 10 minutes in vacuum-assisted steam sterilizers.

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ATCC is a registered trademark of American Type Culture Collection.