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
Quality certification

Bionova® BT220

Esterilización por Vapor / Steam sterilization

Geobacillus stearothermophilus ATCC 7953

LOT

Val

Exp

Población / Population

UFC / CFU

Valor D (121 °C)/ D - value

min.

Tiempo sobrevida / Survival time

min.

Survival time = (log₁₀ labeled population - 2) x labeled D-value

Tiempo de muerte / Kill time

min.

Kill time = (log₁₀ labeled population + 4) x labeled D-value

Valor D (135 °C)/ D - value

min.

Tiempo sobrevida / Survival time

min.

Survival time = (log₁₀ labeled population - 2) x labeled D-value

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Parameters determined at the time of manufacture according to ISO 11138-1: 2006, ISO 11138-3: 2006 e IRAM 37102: 1999 (Parts 1 and 3) standards. The values shown are reproducible only under the same conditions under which they were determined.

ISO and USP Compliant

ATCC is a registered trademark of American Type Culture Collection.

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Composition
Each tube contains a population of *Geobacillus stearothermophilus* ATCC 7953 spores soaked on a filter paper strip (spores carrier). It also has purple colored culture medium contained within the glass ampoule, over the filter next to the carrier, on the base of the tube.

Product description
BT220 Rapid Readout Fluorescence System Biological Indicators for Steam sterilization processes have been designed for quick and easy monitoring of vacuum assisted and gravity displacement steam sterilization cycles at 121-135 °C. They consist of a plastic tube, a special filter, a spores carrier and a glass ampoule with culture medium. The tube has a plastic cap with holes and a barrier permeable to steam.

Rapid readout: 3 hours
The rapid readout must be carried out in the Bionova® IC10/20FR Reader Incubator or in a similar reader incubator that meets the features that are described next. The Reader-incubator must be capable of reading the fluorescence emission of the product, resulting from breaking a specific substrate on the spores carrier. Fluorescence is produced when the reader stimulates the carrier with UV light at 360 nm. Final readout of negative results is readily available after 3 hours of incubation. Fluorescence is a direct measure of the germination and growth of *Geobacillus stearothermophilus* spores which have survived the sterilization process (positive result). Furthermore, a failure in the sterilization process can also become evident by the color change of the culture medium. Due to the high sensitivity of the fluorescence results at 3 hours, conventional incubation for color change of BT220 Rapid Readout Fluorescence System Indicator is not an advantage.

WARNING: Do not use the Bionova® BT220 Biological Indicators to control EO, Dry Heat, Formaldehyde or other sterilization processes. Do not reuse the biological indicators.

Storage
Store at 10-30 °C, 30-80 % Relative Humidity.
Do not freeze.
Do not store biological indicators near sterilizing agents or other chemical products.

Instructions for use
1. Identify the BT220 Bionova® Biological Indicator by writing the sterilizer number (in case of having more than one sterilizer), load number and processing date on the label.
2. Pack the biological indicator along with materials to be sterilized in an appropriate package according to recommended sterilization practices. Place the package in those areas which you consider a priori most inaccessible for the sterilizing agent (Steam). Generally, the center of the load and areas near the door are problematic.
3. Sterilize as usual.
4. After the sterilization process has finished, open the sterilizer door, wait five minutes and remove the biological indicator from the package. **CAUTION:** Wear safety glasses and gloves when removing the Bionova® BT220 Biological Indicator from the sterilized package. **WARNING:** Do not crush or handle the biological indicator excessively, since this might cause the glass ampoule to burst.
5. Let the biological indicator cool down until it reaches room temperature.
6. Check the Chemical Indicator on the label of the biological indicator. A color change to brown indicates that the biological indicator has been exposed to steam at 121 °C or more. **IMPORTANT:** this color change does not evidence the process effectiveness to achieve sterility. If the Chemical Indicator color has not changed, check the sterilization process.
7. Press the lid to seal the tube. Then, break the ampoule contained in the biological indicator. This may be done in 3 different ways:
A- With an ampoule crusher.
B- Manually. To do this, take the tube by placing your index finger and thumb on the space between the cap and the line of chemical indicator and then press the tube.
C- With the ampoule crusher placed within the top of the incubator's incubation area.
Then shake the tube down vigorously, with movements similar to those performed to lower the temperature in a mercury thermometer, until the medium reaches the base of the tube and soaks the spores carrier entirely. Finally, place the biological indicator in the incubator. **IMPORTANT:** Use a non-sterilized biological indicator as a positive control in order to ensure that correct incubation conditions were met; capability of media to promote rapid growth; viability of spores have not been altered due to improper storage temperature; humidity or proximity to chemicals and proper functioning of IC10/20FR Bionova® Incubator. Both the positive control indicator and the processed indicator should belong to the same batch.
8. Incubate the processed biological indicator and the indicator used as positive control for a maximum of 3 hours at 60±2 °C for rapid readout.
Fluorescence detection by the reader (stimulation 340-380 nm / emission 455-465 nm) means a failure in the sterilization process. If no fluorescence is detected at the end of the 3-hour incubation, then, the result is negative (the sterilization process has been effective). The indicator used as positive control must be detected as such by the reader. It is good practice to incubate a positive control for a visual color change. Record the positive results and discard the biological indicators immediately as it is indicated below.
WARNING: Do not reuse the sterilizer until the biological indicator test result is negative.

Visual confirmation: 48 hours
Optionally, you can perform a visual color-change confirmation after a 48-hour incubation. If the sterilization process has not been successful, culture medium will turn to yellow during incubation at 60 °C, indicating the presence of living spores. If sterilization was successful, culture medium will remain purple after the incubation process.
A definitive negative result is obtained after the 48-hour incubation.
The positive control must show a color change from purple to yellow for results to be valid.

Visual readout: 7 days
A 7-day readout is optional and not intended to be routinely performed. This is an initial validation of the 3-hour reading. The fluorescence results are compared to the 7-day visual readings. The sensitivity of the system is determined as the difference between the 7 days positive indicators and false negative indicators (negative fluorescence reading and visually positive) in relation to the positives at 7 days (*). Based on the claimed ≥ 97% readout reliability of the 3-hour biological indicator, there is no advantage in incubating the indicator beyond 3 hours. **NOTE:** If 7-day readout is made, a humidified incubator will be required to prevent medium from drying out.

Disposal
Discard biological indicators after use according to your country's healthcare and safety regulations. The positive biological indicator can be autoclaved at 121 °C for at least 20 minutes, or at 132 °C for 15 minutes in a gravity displacement steam sterilizer, or at 134 °C for 10 minutes in a vacuum assisted steam sterilizer.

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