PCD 20-C

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PCD 20-C

Process Challenge Device.

Test pack for Steam sterilization processes.

Usage

Detects inadequate air removal and steam penetration in dynamic-air-removal (pre-vacuum) and gravity displacement steam autoclaves.

Condition: 132/134 °C in dynamic-air-removal steam sterilization processes ≥ 4 minutes and 121 °C in gravity displacement steam sterilization processes ≥ 30 minutes.

Applicable Regulation

EN ISO 11140-1:2009, ISO 11138-1:2006, ISO 11138-3:2006, ANSI/AAMI ST79

Classification

Class 1, according to risk.

Authorization

Designed under Quality Management Systems ISO 9001:2008 and ISO 13485:2012. ANMAT (Argentinean National Administration of Drugs, Food and Medical Technology) PM 1614-4.

Characteristics

Porous cards system holding a Self-contained Biological Indicator and IT26-C moving-front Class 5 Integrator Indicator. The system also contains a self-adhesive Record Card, PCDBI-C-RC, where data on sterilization cycle may be written.

Cardboard box (10.5 cm x 6.6 cm x 1.8 cm) with a self-adhesive Class 1 Process Indicator.

Process Indicator:

Initial color: **blue**

Final color: **dark grey or black**

IT26-C Integrator Indicator (within the challenge device):

The migration is visible through a zone marked ACCEPT or REJECT, thus indicating whether sterilization conditions were met.

Integration condition is calibrated against the kill time of a 10^6 *G. stearothermophilus* ATCC 7953 spores' population, calculated in BIER (Biological Indicator Evaluator). Conditions: saturated steam at 121 °C, 128 °C, 135 °C.

Self-contained Biological Indicator (within the challenge device):

$\geq 10^5$ *G. stearothermophilus* ATCC 7953 spores per vial in paper carrier. Glass ampoule with growth indicator medium. Label with chemical indicator line printed with Steam reactive ink (Color change: pink to brown).

Final results: incubation for 24 hours at 60 °C (maximum).

NOTE: If sterilization process was not successful, indicator medium will turn yellow after incubation, thus indicating the presence of living *G. stearothermophilus* spores. If sterilization process was successful, indicator medium will remain purple after incubation.

Environmental conditions during manufacture

Room temperature 15-30 °C, RH 30-80 %.

Storage conditions

Best stored in original box under the next conditions: T = 10-30 °C, RH 30-80 %.

Transportation conditions

Storage conditions should be strictly followed. Carry in closed and reinforced boxes in order to avoid damages. Transportation of this product does not represent a risk for health.

Shelf-life

24 months.

Packing

Packing information: product code and description, operating instructions, classification according to regulation, process for intended use, batch number, expiration date and manufacturer information.

Labelling

Information only on product's box.

Possible target markets

Healthcare, Food, Pharmaceutical and Medical Industries.

Other important information (optional)

Read product's directions for use thoroughly before use.

Note: When necessary, limits for each product's parameters should be specified.

Description of identified and/or predictable risks which could represent a risk for the product

Do not store the product near sterilizing agents.

Do not expose this product to EO, Dry Heat, Radiation or any sterilization process other than Steam.

The positive biological indicator can be autoclaved at 121 °C for at least 20 minutes, at 132 °C for at least 15 minutes in gravity displacement steam sterilization, or at 134 °C for at least 10 minutes in vacuum steam sterilizer.