
Indicadores Biológicos Biological Indicators

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Survival time = not less than (log₁₀ labeled population - 2) x labeled D-value

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Kill time = not more than (log₁₀ labeled population + 4) x labeled D-value

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Parameters determined at time of manufacture according to ISO 11138 (Parts 1) and IRAM 37102 (Parts 1) standards. The shown values are reproducible only under the same conditions under which they were determined.

Lic. Adrián J. Rovetto
Director Técnico
Quality Assurance Director

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BT80

Biological Indicators For Moist Heat decontamination

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Composition

Each tube contains a population of spores of *Bacillus atrophaeus* ATCC 9372 soaked on a carrier. It also has a growth indicator medium contained in the glass ampule.

Product description

Bionova® BT80 biological indicator is specifically designed for the monitoring of moist heat decontamination processes.

If the decontamination process was not successful, the indicator medium will turn to yellow after incubation at 37±2 °C, thus indicating the presence of living *Bacillus atrophaeus* spores.

If the decontamination process was successful, the indicator medium will remain the original color after incubation. The final readout should be carried out after 48 hours of incubation at 37±2 °C.

Warning!

Do not use BT80 biological indicator for monitoring steam sterilization cycles, dry heat, chemical vapor, ethylene oxide, radiation or other different sterilization process for which it was designed. Do not reuse biological indicators.

Storage

Store it in a dark place at temperatures between 10 - 30 °C, 30 - 80 % relative humidity.

Do not freeze.

Do not store these biological indicators near sterilizing agents or other chemical products.

Directions for use

1. Identify the Bionova® BT80 indicator tube by writing the machine number (in case there is more than one), load number, and processing date on the indicator label.

2. Place the biological indicator along with the items to be decontaminated in an appropriated package according to recommended decontamination practices. Place this package in those areas which a priori you consider most inaccessible for decontamination agent (Moist Heat). Generally a problematic area is the center of the load.

3. Proceed as usual.

4. After the decontamination process has finished, allow the load to cool until it is at/below room temperature and remove the biological indicator from the load for its incubation.

5. Crush the glass ampoule contained in the biological indicator and incubate at 37±2 °C.

IMPORTANT: Use a non-processed biological indicator as positive control each time a processed indicator is incubated. The positive control ensure that correct incubation conditions were met.

6. Incubate the processed biological indicator and the indicator used as positive control for a maximum of 48 hours at 37±2 °C. Reading should be performed at convenient intervals of 10 hours.

A color change to yellow of the growth indicator medium means a decontamination process failure has occurred. If after 48 hours there is no color change in the processed indicators, a final negative result is obtained, thus the decontamination process was successful. The positive control indicator should show a yellow color change for the results to be valid.

Record the positive ones and discard them immediately as is shown below.

WARNING! Do not reuse the decontamination machine until the biological indicator test results are negative (process indicator remain the original color).

Disposal

Discard biological indicators after use according to your country's healthcare and safety regulations. The positive biological indicator can be autoclaved at 121 °C for at least 20 minutes, or at 132 °C for 15 minutes in a gravity displacement steam sterilizer, or at 134 °C for 10 minutes in a vacuum assisted steam sterilizer.

TERRAGENE

Industria Argentina - Made in Argentina

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Fabricado por: Terragene S.A. - Güemes 2879 - (2000) Rosario - Santa Fe - Argentina