
BT90 Biological Indicator

Rev. 13 / 07.2017

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Biological Indicator

For Plasma or Vaporized Hydrogen Peroxide sterilization

Composition

Each tube contains a population of *Geobacillus stearothermophilus* ATCC 7953 spores laid on the inner surface of the plastic tube base. It also has a purple growth indicator medium contained in the glass ampoule, on a plastic filter.

Description
BT90 Rapid Readout Colorimetric System Biological Indicators have been designed for quick and easy monitoring of Plasma or Vaporized Hydrogen Peroxide sterilization processes. The system contains a tube with a conic base and a special filter that keeps the glass at the top of the tube. An easy readout can be carried out using the reference color chart. The tube has a plastic cap with holes and a barrier permeable to Hydrogen Peroxide. If the sterilization process was not successful, the indicator medium will change from purple to greyish-green and finally to yellow after incubation at 60±2 °C, thus indicating the presence of living *Geobacillus stearothermophilus* ATCC 7953 spores.
If the sterilization process was successful, the indicator medium will remain purple after incubation. The final readout should be performed after 8 hours of incubation at 60±2 °C.

WARNING: Do not use BT90 Biological Indicators to control EO, Steam, Dry heat, Formaldehyde or other sterilization process. Do not re-use biological indicators.

Storage
Store in a dark place and at temperatures between 10-30 °C, 30-80 % Relative Humidity.
Do not freeze.
Do not store these biological indicators near sterilizing agents or other chemical products.

Directions for use
1. Identify the Bionova® BT90 indicator by writing the sterilizer number (in case of having more than one), load number, and processing date on the indicator label.
2. Pack the biological indicator along with the materials to be sterilized in an appropriate package according to recommended sterilization practices. Place this package in those areas which a priori you consider most inaccessible for the sterilizing agent. Generally a problematic area is the center of the load and near the door.
3. Sterilize as usual.
4. After the sterilization process has finished, open the sterilizer door and remove the biological indicator from the package. **WARNING!** Do not crush or handle the biological indicator excessively, since this might cause the glass ampoule to burst.
5. Check the chemical indicator on the label of the biological indicator. A color change to green confirms that the biological indicator has been exposed to Vaporized Hydrogen Peroxide. **IMPORTANT:** This color change does not indicate that the process was sufficient to achieve sterility. If the chemical indicator is unchanged, check the sterilization process.
6. Break the ampoule contained in the biological indicator either with an ampoule crusher or manually. If you use an ampoule crusher, you must first press the lid down to seal the tube before using the device (fig. A). If you do it manually, you must hold the tube by the space between the cap and the chemical indicator line, press to break the ampoule and then press the lid down to seal the tube (fig. B).
7. Incubate the processed biological indicator as well as the control used as an ampoule crusher or manually, with movements similar to those performed to lower the temperature in a mercury thermometer, in order to collect the liquid in the conical base of the tube (fig. C). **IMPORTANT:** Use a non-sterilized biological indicator as a positive control to ensure that correct incubation conditions were met.
8. Incubate the processed biological indicator as well as the positive control for a maximum of 8 hours at 60 ± 2 °C. Color change of the growth indicating medium from purple to greyish green, indicates failure in the sterilization process. If after 8 hours, color change of the processed indicators is not observed, the result is negative (the sterilization process was satisfactory). Color of the indicator used as positive control should turn from purple to greyish green and finally to yellow for results to be considered valid.

DETECTION METHOD. THE SYSTEM ALLOWS AN EARLY DETECTION OF THE POSITIVE BIOLOGICAL INDICATORS BY A GRADUAL COLOR CHANGE OF THE CULTURE MEDIUM.
Readout should be carried out by using the reference color chart, under white light. Yellow light might complicate the readout. The chart enables a simple detection of the positive biological indicators, allowing a final readout during the first 8 hours of incubation with a sensitivity of 97%. For an easy readout, shake the IB down to collect the liquid at the bottom of the tube. Record the results and immediately discard as indicated below. To remove the label of the tube for data registration, the tube cap should be lifted up against the original position and the label should be peeled off. If the cap is not 3 mm high, the label could be damaged.

NOTE: BT90 readout/incubation time is 8 hours at 60±2 °C. It is not recommended to extend the incubation time beyond 8 hours, since it is likely that sterilizing agent residues could affect the color of the culture medium, bringing unreliable results. However, if you wish to extend the incubation time beyond 8 hours, it is advisable to use a humidified environment in order to prevent complete evaporation of the culture medium of the indicator. **WARNING!** Do not use the sterilizer again until biological indicators results are negative.

Disposal
Discard biological indicators after use according to your country's healthcare and safety regulations. The positive biological indicators can be autoclaved at 121 °C for at least 20 minutes, or at 132 °C for 15 minutes in a gravity air displacement steam sterilizer, or at 134 °C for 10 minutes in a vacuum assisted steam sterilizer.

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ATCC is a registered trademark of American Type Culture Collection