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# Cleaning Challenge Device for flexible endoscopes

For Automated Endoscope Reprocessors

EN

Endoscopes are very complex reusable medical devices. They consist of several channels and tubes which can be 2 meters long, with variable internal diameters. Endoscopes get in contact with human mucous membranes and tissues, consequently, high-level cleaning and disinfection is required after its use, in order to avoid pathogen transmission. Automated Endoscope Reprocessors (AERs) are the preferred method for endoscope washing and disinfection given that they provide a standardized method with minimal human intervention, together with the fact that washing and disinfection cycles are highly reproducible. An optimal result depends on the adequate combination of water pressure and quality, chemical agents' concentration, temperature and cycle duration. Monitoring the efficacy of cleaning processes is crucial in order to ensure compliance with the standards of cleanliness.

## Product description

Chemdy® CCDER1 and CCDER2 devices consist of a perforated plastic capsule able to hold the cleaning indicator CDWE which contains a mixture of specially combined colored components. The capsule is connected to a high resistance plastic tube of 2 meters long and 1 mm (CCDER1) or 2 mm (CCDER2) internal diameter. As well, the tube has a luer-lock connector, which enables the connection to Automatic Endoscope Reprocessors terminals. One luer-lock male and female terminal is included. Select the most appropriated terminal for your AER.

## Instructions for use

- 1-Place a Chemdy® CDWE indicator strip inside the capsule, with the indicator ink towards the perforated side.
- 2-Screw the capsule lid.
- 3-Connect the Chemdy® CCDER device to a proper AER terminal.
- 4-Connect the endoscope terminals to proper AER terminals.
- 5-Run the selected washing program.
- 6-Once the process is finished, disconnect the Chemdy® CCDER device from the AER.
- 7-Remove the indicator strip and evaluate the result.
- 8-Register the adhesive strip on your record book.

## Results interpretation

Complete elimination of the washing ink indicates that cleaning efficacy was adequate. If ink residues are observed, the endoscope cannot be considered clean. Failures during a cleaning cycle can be a consequence of inadequate chemical agents' concentration, insufficient water pressure, poor water quality, temperature or cycle duration. All these events diminish cleaning efficacy and it is not recommended to use the endoscope until a satisfactory result is obtained.

**NOTE:** Endoscope reprocessing monitoring can be complemented with chemical indicators for peracetic acid as disinfectant, the Integron® IT400 and IT401 series. Chemdy® CDWE and CCDER devices allow monitoring endoscope cleaning efficacy but not disinfection level. For disinfection level monitoring use the biological indicator Bionova® BT400, following the Instructions for use of this product.

## Additional Technical Information

In each cycle, the user should check that the connectors and the plastic tube are in perfect conditions. In case that the plastic tube is blocked or has leakages, the complete device should be replaced. If the connector is damaged and its replacement is justifiable, contact a representative and request a replacement. A luer-lock connector is attached.

**CAUTION:** Terragene® does not guarantee the proper functioning of Chemdy® CCDER devices if used with washing indicators of other manufacturers.

Results obtained with the device Chemdy® CCDER are only valid for the wash cycle under analysis. Conclusions for previous or subsequent cycles should not be drawn.

Terragene® does not guarantee the proper functioning of Chemdy® CCDER if used for different washing processes than the ones it was designed for.

## Storage

Store the Chemdy® CCDER devices in a dry place, protected from light, at a temperature from 10-30 °C and 30-80% Relative Humidity.

**NOTE:** incorrect storing of the product may affect the stability thereof and the results obtained.

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