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# Multiple Cleaning Challenge Devices for flexible endoscopes

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PT

## For Automated Endoscope Reprocessors

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Endoscopes are very complex reusable medical devices. They consist of several channels and tubes which can be 2 meters long, with variable internal diameters. Endoscopes get in contact with human mucous membranes and tissues, consequently, high-level cleaning and disinfection is required after its use, in order to avoid pathogen transmission. Automated Endoscope Reprocessors (AERs) are the preferred method for endoscope washing and disinfection given that they provide a standardized method with minimal human intervention, together with the fact that washing and disinfection cycles are highly reproducible. An optimal result depends on the adequate combination of water pressure and quality, chemical agents' concentration, temperature and cycle duration. Monitoring the efficacy of cleaning processes is crucial in order to ensure compliance with the standards of cleanliness.

### Product description

Chemdye® CCDER122 device consist of three perforated plastic capsules which are able to hold the cleaning indicator CDWE and are connected to 2 meters long high resistance plastic tubes. One of the tubes has 1 mm of internal diameter and the other two have 2 mm of internal diameter. As well, each tube have a luer-lock connector which enable the connection to Automatic Endoscope Reprocessor terminals. One luer-lock male and female terminal is included for each tube. Select the most appropriated terminal for your AER.

### Instructions for use

- 1-Place a Chemdye® CDWE indicator strip inside each capsule's lid, with the indicator ink towards the perforated side.
- 2-Screw the capsule's lid.
- 3-Connect the Chemdye® CCDER122 device to proper AER terminals.
- 4-Connect the endoscope terminals to proper AER terminals.
- 5-Run the selected washing program.
- 6-Once the process is finished, disconnect the Chemdye® CCDER122 device from the AER.
- 7-Remove the indicator strip and evaluate the result.
- 8-Register the adhesive strip on your record book.

### Results interpretation

Complete elimination of the washing ink indicates that cleaning efficacy was adequate. If ink residues are observed, the endoscope cannot be considered clean. Failures during a cleaning cycle can be a consequence of inadequate chemical agents' concentration, insufficient water pressure, poor water quality, temperature or cycle duration. All these events diminish cleaning efficacy and it is not recommended to use the endoscope until a satisfactory result is obtained.

**NOTE:** Endoscope reprocessing monitoring can be complemented with chemical indicators for peracetic acid as disinfectant, the Integron® IT400 and IT401 series. Chemdye® CDWE and CCDER devices allow monitoring endoscope cleaning efficacy but not disinfection level. For disinfection level monitoring use the biological indicator Bionova® BT400, following the Instructions for use of this product.

### Additional Technical Information

In each cycle, the user should check that the connectors and the plastic tube are in perfect conditions. In case that the plastic tube is blocked or has leakages, the complete device should be replaced. If the connector is damaged and its replacement is justifiable, contact a representative and request a replacement. Three luer-lock connectors are attached.

**CAUTION:** Terragene® does not guarantee the proper functioning of Chemdye® CCDER122 devices if used with washing indicators of other manufacturers.

Results obtained with the Chemdye® CCDER122 device are only valid for the wash cycle under analysis. Conclusions for previous or subsequent cycles should not be drawn.

Terragene® does not guarantee the proper functioning of Chemdye® CCDER122 if used for different washing processes than the ones it was designed for.

### Storage

Store the Chemdye® CCDER122 devices in a dry place, protected from light, at a temperature from 10-30 °C and 30-80% Relative Humidity.

**NOTE:** incorrect storing of the product may affect the stability thereof and the results obtained.

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