Cleaning Performance Test

BASKET ZONE: water pressure, spray system performance, detergent chemistry and concentration, temperature.

COVERED ZONE: detergent chemistry and concentration, water quality, temperature, cycle duration, etc.

Results reference guide (✓OK / ✗ FAIL)

Visual results	BASKET ZONE influence results			COVERED ZONE influence results			
	Location inside washer	Water pressure	Mechanical functions (detergent injection pump, spray arms, etc.)	Temperature	Exposure time	Detergent (type, dose)	Water quality
Processed and correct	✓	✓	✓	✓	✓	✓	✓
✗	✓	✓	✓	✗	✗	✗	✗
✗	✓	✓	✓	✗	✗	✗	✗
✗	✓	✓	✓	✗	✗	✗	✗
✗	✓	✓	✓	✗	✗	✗	✗
✗	✓	✓	✓	✗	✗	✗	✗
✗	✗	✗	✗	✗	✗	✗	✗
✗	✗	✗	✗	✗	✗	✗	✗
✗	✗	✗	✗	✗	✗	✗	✗
Unprocessed	—	—	—	—	—	—	—

NOTE: These results could be a consequence of failure in one or a combination of variables. When detecting the failed variable/s, make the corrective actions. If the problem persists, contact the technical service of your washing machine. The higher the quantity of ✗, the higher the influence of that variable on the failure of the cycle. Reference colors exhibited in printed prospectus as well as those displayed in our website and soft copies of other documents, do not necessarily represent the real color as shown in the actual indicators.

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ES

Cleaning Indicators
For cleaning performance tests

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PT

Control of the process of cleaning and washing of medical instruments and materials is highly important since the outcome affects subsequent disinfection and/or sterilization processes success. Cleaning control is particularly relevant in automatic washing machines, where even though process parameters may be well controlled, the cleaning final result depends largely on uncontrolled factors. Washing efficacy is influenced by detergent, water temperature, time and washing penetration. A good cleaning result is achieved by an appropriate combination of these variables.

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Indicators Classification

Code	Sample color	Challenge level	Annex of standard	Standard reference test soil content
CDWA3	Blue	High	G	Semolina pudding
CDWA3	Azul	Alto	G	Budin de sémola
CDWA4	Fuchsia	Very high	Q	Sheep blood, water soluble wallpaper paste, egg, black ink

NOTE: It is recommendable to test the 2 types of CDWA indicators along with the instrument, in order to select the one which is more appropriate for the washing machine performance. It is advisable to start the washing control validation process using the most rigorous CDWA product. The indicator that washes completely during the washing/dishwashing program could be used for routine monitoring. It should be noticed that the choice of indicator will depend on the hardness of dirt adhered to the object to be washed. Therefore, the CDWA indicator selected for routine monitoring of a particular washing machine could be changed according to the hardness of the soil on the load.

Once the CDWA indicator appropriate for the washing requirements has been selected, we recommend placing several identical indicators at different positions of the washing machine, so as to monitor the performance of all areas of the machine during the washing process.

Characteristics

CDWA indicators have a space for recording information and are placed within a special holder (CDWA holder, see figure) that can be fixed on a tray. Holder design mimics a realistic challenge of thoroughly cleaning hard-to-reach areas to the washing/dishwashing process. When CDWA is placed on the holder, the top half of its surface is obscured (COVERED ZONE, see the figure), thus simulating the effect of a surgical instrument with water stream hitting it indirectly. The other half of the test soil is covered by a mesh (BASKET ZONE in the figure), thus posing a difficulty for water stream to reach it, mimicking shadowed or occluded surfaces.

CDWA indicators provide constant cleaning characteristics compared to the use of natural samples based assays which might suffer modifications of adhesion and polarity properties during storage time.

The test system provides high confidence and high reliability to reveal defects in the cleaning process.

CDWA indicator could also be used for routine monitoring of ultrasonic washers, which operate at frequencies equal or greater than 35 kHz, between 20°C and 65°C (validation on each washer is required).

WARNING: The indicator should be placed in the middle of the water bath, otherwise it could cause erroneous results.

Instructions for use

- Choose the CDWA indicator according to the desired difficulty.
- Place the indicator on the CDWA holder.
- Place the system on the tray along with the material.
- Once the cleaning cycle is over, remove the indicator.
- Analyze the indicators CDWA with the table of reference and register the results. Ante un resultado desfavorable, realice las acciones correctivas pertinentes y repita el proceso de lavado.

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NOTE: For the holder model in which the test soil is completely covered by a mesh, the effect of washing/dishwashing process is homogeneous, unlike to the latest holder design with mesh only in the top half, which allows simulating two different effects. Therefore, if the full mesh holder is used, it should be noted that the result shown in the reference table for the top half of the test soil will be repeated in the bottom half.

Observations

Technical specifications outlined in ISO/TS 15883-5:2005 standards describe different types of test samples having different properties and hardness facing cleaning, but no recommendation is made about which sample type should be used. Currently, none of the samples is suggested as a reference or standard because no method has been defined to compare the properties of the materials used in these tests. For this reason, Terragene® provides two washing indicators with different levels of response and adhesion.

Precautions

- Do not reuse CDWA indicators.
- Do not touch the indicator ink with your hands.

Storage

CDWA indicators should be stored at temperatures between 10-30 °C, protected from light and moisture, within its original packaging.

NOTE: in case of remaining traces of indicator colored spot after the washing process, do not store the used indicator under direct light since this may deteriorate its coloring.

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CDWA cleaning system

CDWA cleaning system offers two different levels of cleaning with different adhesion characteristics, thus allowing different challenge levels to test cleaning efficacy. The indicators efficiency is validated by comparing against the test soils detailed in ISO/TS 15883-5:2005 standard annexes.

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