IT27-3.5S - IT27-3.5SAD - IT27-3.5L - IT27-5.3S
IT27-5.3SAD - IT27-5.3L - IT27-7S - IT27-7SAD
IT27-9S - IT27-18S

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Results reference guide

Unprocessed

Failures
One or more parameters of sterilization were not achieved. The materials cannot be considered sterile and must be processed again. These results are usually obtained when steam quality is poor or adequate vacuum has not been generated

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Type 6 Indicator according to ISO 11140-1:2014

Product description

The Integron® Emulator Indicators are designed to react in steam sterilization processes inside of the loads, ensuring an adequate control of the effectiveness of sterilization processes (steam, time, quality of steam). The purple ink was developed to change to green when the process reaches the variables set out in the emulator.

Characteristics

IT27-3.5S, IT27-3.5SAD, IT27-5.3S, IT27-5.3SAD, IT27-7S, IT27-7SAD, IT27-9S and IT27-18S emulator indicators consist of a paper strip of 72 mm x 23 mm. IT27-3.5L and IT27-5.3L emulator indicators consist of a paper strip of 142 mm x 23 mm.

IT27-3.5SAD, IT27-5.3SAD and IT27-7SAD emulator indicators are printed on a self-adhesive support which allows their easy record.

Indicators 100 % metals free.

Indications for use

Use Integron® Emulator Indicators for monitoring Steam sterilization processes.

Emulation condition:

IT27-3.5S	3.5 min.	134 °C	
IT27-3.5SAD	3.5 min.	134 °C	
IT27-3.5L	3.5 min.	134 °C	
IT27-5.3S	5.3 min.	134 °C	15 min. 121 °C
IT27-5.3SAD	5.3 min.	134 °C	15 min. 121 °C
IT27-5.3L	5.3 min.	134 °C	15 min. 121 °C
IT27-7S	7 min.	134 °C	20 min. 121 °C
IT27-7SAD	7 min.	134 °C	20 min. 121 °C
IT27-9S	9 min.	134 °C	
IT27-18S	18 min.	134 °C	

Precautions

Do not expose these products to EO, Dry Heat, or any sterilization processes other than Steam. Do not reuse the sterilizer until the emulator's color reaches the emulation condition.

Instructions of use

- Place an Integron® Emulator Indicator within each package, bag or tray that will be sterilized by Steam.
- In rigid containers place an emulator indicator in every corner, or at least two diagonally or opposite corners, to ensure that the steam has reached the most inaccessible areas.
- After completing the process, analyze the results. The indicator must change to the reference color for considering that emulation conditions were reached.

Storage and shelf life

- Indicators should be stored in dry places, protect from light at temperatures between 10-30 °C, 30-80 % Relative Humidity. Do not wet. Do not storage close to sterilizing agents.
- Indicators have an expiration date of 5 years from the date of manufacture when stored in the recommended conditions. Both expiration date and batch number are specified on package label and on each unit of product. Do not use indicators after their expiration date.
- After use, the indicator will not change in 1 year if stored in temperature and relative humidity previously indicated.

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