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ES

## Thermodisinfection and cleaning indicators

For moist heat disinfection processes in washer-disinfectors

## Product description

After cleaning, surgical instruments must be disinfected to ensure safe subsequent handling and processing. Thermal disinfection with moist heat, is the most common method for disinfection of medical devices in the hospital setting. Disinfection effects can be defined by the  $A_{v0}$  value, a physical parameter specified in ISO 15883-1, which denotes the inactivation of microorganisms. A particular  $A_{v0}$  value can be achieved with the most diverse temperature/time combinations. For processing surgical instruments, an  $A_{v0}$  value of 600 is specified in ISO 15883-2, corresponding to a hold time of 1 minute at 90 °C, ensuring lethality of bacteria, mycobacteria, fungi and heat-sensitive viruses. If efficacy is also to be ensured against heat-resistant viruses, e.g. hepatitis B, many authorities in the field of disinfection, like the Robert Koch Institute, recommend an  $A_{v0}$  value of 3000, corresponding to a temperature of 90 °C with a hold time of 5 min. A higher safety approach from this institute comprises a cycle at 93°C, 10 minutes, a program known as the 'BGA/RKI' or epidemic program. Integron® Indicators have been designed to react in moist heat disinfection processes in washer-disinfectors, providing an accurate, convenient method of routine control, and assuring disinfection efficiency. The green indicating ink was developed to turn to purple when temperature and time conditions are met.

## Characteristics

Integron® Indicators consist of one-sided-laminated strips of 30 mm x 70 mm printed with green reactive ink, designed with state-of-the-art technology to withstand conditions present during moist heat disinfection processes.

## 100% Toxic Heavy Metals free

## Indications for use

Use the correspondent Integron® Chemical Indicator for monitoring thermal disinfection processes in washer-disinfectors according to the following programs:

IT27W-1: 90°C, 1 minute.  $A_{v0}=600$   
IT27W-5: 90°C, 5 minutes.  $A_{v0}=3000$   
IT27W-10: 93°C, 10 minutes. BGA/RKI/Epidemic Program

## Precautions

Do not expose this product to Steam, EO, Dry Heat, or any sterilization or disinfection processes other than Moist Heat Disinfection Processes in washer-disinfectors. Do not reuse.

Warning! The results are not reproducible when immersing the indicators in a thermal bath.

## Instructions for use

- 1- Place the indicator on the holder.
  - 2- Place the system on the tray along with the material.
  - 3- Run the selected disinfection program and once finished remove the indicator.
  - 4- Visually analyze indicators and record results. The indicator must change to the reference color for considering that the desired disinfection conditions were met. If unsatisfactory results are obtained, perform the actions indicated therein and repeat the disinfection process.
- Use one Integron® indicator for each cycle.

## Interpretation of results

The color change of the indicating ink from green to purple indicates that the thermal disinfection conditions were attained, thus allowing a quick and easy assessment of the success of the process. If no color change is observed, the load cannot be considered disinfected. Failures during a disinfection cycle can be a consequence of inadequate low temperature or cycle duration, which lowers its efficiency. It is not recommended to release the load until a satisfactory result is obtained.

## Storage and shelf life

Indicators should be stored in dry places, away from light, at temperatures between 0°C and 4°C, 30-80% relative humidity. Do not wet. Do not store close to sterilizing agents.

Indicators have an expiration date of 5 years from the date of manufacture when stored in the recommended conditions. Both expiration date as well as the batch of the product are specified on the package label and on each unit of product. Do not use indicators after their expiration date.

EN