## PCD-A-P1

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## PCD-A-P1

For Hydrogen Peroxide sterilization

EN

Chemical Strips (Type 2 Chemical Indicator according to ISO 11140) for use in Helix-PCD for Hydrogen Peroxide.

## Instructions for use

1-A. Remove the threaded cap from the metal or plastic capsule. Only for metal capsule Helix devices, remove the bolt attached to the cap (the bolt should NOT be used in Hydrogen Peroxide processes).

1-B. Place the chemical indicator in the lumen of the capsule.

**IMPORTANT:** Avoid folding the chemical indicator in the sectors with chemical inks.

2. Place the lid of the capsule with moderate pressure. **ATTENTION!** Do not place the system inside the cloth bag, since it can interfere with the process. The cloth bag must be used in steam processes.

**NOTE:** Terragene® suggests using the system Chemdye® Helix-PCD designed for small sterilizers (CDPCD2X025).

The system Chemdye® Helix-PCD should be placed with the camera completely empty.

3. For Hydrogen Peroxide processes, start the selected program and test the process.

Initial color: purple

Final color: green

4. After the period of ventilation, remove the indicator and evaluate the result. If the 4 bars changed to the final color, the result is successful, the sterilizer is working correctly. If one bar is different to the final color, it shows the presence of residual air in the sterilizer, therefore an equipment malfunction.

We suggest retrying the test once again, before sending the equipment to technical service.

Color change of Hydrogen Peroxide indicator strip:

 Initial condition

 Final condition

## Storage

Indicators should be stored in dry places, protect from light at temperatures between 10-30 °C, 30-80 % Relative Humidity.

Do not wet.

Do not store close to sterilizing agents.

Do not use indicators after their expiration date.

Do not reuse.

**NOTE:** Reference colors exhibited in printed prospects as well as those displayed in our website and soft copies of other documents, do not necessarily represent the real color as shown in the actual indicators.

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