# BTD40-3/6

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# BTD40-3/6

## Spore Discs Biological Indicator

For Ethylene Oxide (EO) and Dry Heat sterilization processes.

### Usage

Monitoring EO and Dry Heat Sterilization Processes. For industrial use only.

### Applicable Regulation

ISO 11138-1:2006, ISO 11138-2:2006 and ISO 11138-4:2006; IRAM 37102-1:1999 and IRAM 37102-2:1999.

### Classification

Class 1, according to risk.

### Authorization

Designed under Quality Management System standards ISO 13485:2003/NS-EN ISO 13485:2012. ANMAT (Argentinean National Administration of Drugs, Food and Medical Technology) PM 1614-1.

### Characteristics

1/3" (8.46 mm) diameter filter paper disc.

10<sup>6</sup> *Bacillus atrophaeus* ATCC 9372 spores per disc.

Upon completion of the sterilization cycle, the spore disc should be tested in a suitable culture medium at 37°C. We suggest the use of BIONOVA® MC1030 or MC1020 culture media to final reading in 48 hours. For conventional culture medium like TSB, incubation for 7 days at 37°C is recommended.

D-Value<sub>DRY</sub>: Not lower than 2.5 min at 160 °C. Another D-Value is estimated at 170 °C for calculation of Z-Value.

Z-Value: Not lower than 20 °C.

D-Value<sub>EO</sub>: Not lower than 2.5 min, conditions: 600 mg/l EO, 60 % RH, 54 °C for EO.

### Environmental conditions during manufacture

Room temperature 15-30 °C, RH 30-80 %. Sterility conditions are necessary only during the inoculation process performed in laminar flow.

### Storage conditions

T = 10-30 °C, RH 30-80 %, keep in a dark place.

### Transportation conditions

Storage conditions should be strictly followed.

Products should be transported in closed and reinforced boxes in order to avoid damages. Product transportation does not represent any risk for human health.

### Shelf-life

24 months.

**Packing**

100 units per envelope.

Packing information: operating instructions, classification according to regulation, manufacturer information and data on envelope 's label.

Weight per envelope: 14,0 g.

**Labelling**

On product's envelope: product code and description, process for intended use, operating instructions, presentation, batch number, expiration date and manufacturer information.

**Possible target markets**

Industry.

**Other important information (optional)**

It is advisable to incubate at 37 °C in BIONOVA® IC10/20 incubator.

Read product's directions for use thoroughly before use.

NOTE: When necessary, limits of each product's parameters should be specified

**Description of identified and/or predictable risks which could represent a risk for the product**

Do not store the product near sterilizing agents.

Do not expose this product to radiation, H<sub>2</sub>O<sub>2</sub> or any sterilization process other than EO and/or Dry Heat.

The positive biological indicator can be autoclaved at 121 °C for at least 20 minutes, at 132 °C for at least 15 minutes in gravity displacement steam sterilization, or at 134 °C for at least 10 minutes in vacuum steam sterilizer.

Note: Spore discs can be incubated BIONOVA® MC1030 or MC1020 culture media or other culture medium made up of soybean casein and yeast extract.