# CDPCD2X025/P

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# CDPCD2X025/P

## Helix Process Challenge Device

*To evaluate air removal and penetration of sterilizing agent in hollow loads.*

### Usage

To evaluate air removal and penetration of the sterilizing agent (steam / ethylene oxide / formaldehyde /  $\text{VH}_2\text{O}_2$ ) in hollow loads.

Conditions:

STEAM: 3.5 min., 134 °C // 5.3 min., 134 °C // 7.0 min., 134 °C // 9.0 min., 134 °C // 18 min., 134 °C.

OE: 240 min., 54 °C, EO 600 mg/l, 60 % RH.

FORM: LTSF processes.

$\text{VH}_2\text{O}_2$ : plasma or vaporized Hydrogen Peroxide processes.

### Applicable regulation

Less demanding requirements than EN 867-5.

### Authorization

Designed under Quality Management System standards ISO 9001:2008 and ISO 13485:2012.

### Characteristics

Plastic (polypropylene) capsule of great resistance and durability, to hold a chemical indicator designed for this system. The capsule is connected to a high resistance plastic tube of 0.25 m length and 2 mm internal diameter. Special silicone connector coated by external connector.

NOTE: Terragene recommends Helix-PCD of 0.25, 0.5, 0.75 or 1.2 m length for  $\text{VH}_2\text{O}_2$  processes.

Chemical indicators designed for Helix-PCD system: PCD-A-3.5 / PCD-A-5.3 / PCD-A-7.0 / PCD-A-9.0 / PCD-A-18 (STEAM); PCD-A-F1 (FORM); PCD-A-E1 (EO); PCD-A-P1 ( $\text{VH}_2\text{O}_2$ ).

NOTE: If 4 bars change to the final color, the result is successful, the sterilizer works properly. If one or more bars have not turned to the final color, this shows a failure in sterilizing agent penetration, therefore an equipment malfunction.

### Environmental conditions during manufacture

15-30 °C (room temperature), RH 30-80 %.

### Storage conditions

T = 10-30 °C, RH 30-80 %, keep out of direct light.

### Transportation conditions

Storage conditions should be strictly followed. Transport in closed and reinforced boxes in order to avoid damages. Transportation of this product does not represent a risk for health.

### Shelf-life

Helix-PCD: 250 cycles

**Packing**

1 Helix-PCD + 1 Cotton bag

Packing information: product code and description, manufacturer information and data on label.

**Labelling**

Only on product's box: product code and batch number.

**Possible target markets**

Healthcare, Food, Pharmaceutical and Medical industries.

**Other important information (optional)**

Read product's directions for use thoroughly before use.

**Note:** When necessary, limits of each product's parameters should be specified.

**Description of identified and/or predictable risks which could represent a risk for the product**

Do not place the Helix-PCD device within the cotton bag when sterilizing with  $\text{VH}_2\text{O}_2$  and Formaldehyde; since it can interfere with these processes.

Do not store the product near sterilizing agents.

Do not expose this product to chemical agents.