TERRAGENE

Product Name

CDPCD2X45/M – CHEMDYE® TERRAGENE® Helix Process Challenge Device.

Usage

To evaluate air removal and penetration of the sterilizing agent (steam / ethylene oxide / formaldehyde / H₂O₂).

Conditions:

STEAM: 3.5 min., 134 °C // 5.3 min., 134 °C // 7.0 min., 134 °C // 9.0 min., 134 °C // 18 min., 134 °C.

OE: 240 min., 54 °C, EO 600 mg/l, 60 % RH.

FORM: LTSF processes.

H₂O₂: plasma or vapor processes.

Applicable Regulations

EN 285 / EN 867-5

Product Classification

Class 1, according to risk.

Authorization

Designed under Quality Systems ISO 9001:2008 and ISO 13485:2012.

Product characteristics

Metal capsule of great resistance and durability, to hold a chemical indicator designed for this system. The capsule is connected to a high resistance plastic tube of 4.5 m length and 2 mm internal diameter. Two special silicone connectors coated by external connector.

NOTE:

Terragene recommends Helix-PCD of 0.25 m length for H₂O₂ processes.

Chemical indicators designed for Helix-PCD system: PCD-A-3.5 / PCD-A-5.3 / PCD-A-7.0 / PCD-A-9.0 / PCD-A-18 (STEAM); PCD-A-F1 (FORM); PCD-A-E1 (EO); PCD-A-P1 (H₂O₂).

NOTE:

If 4 bars change to the final color, the result is successful, the sterilizer works properly. If one or more bars have not turned to the final color, this shows a failure in sterilizing agent penetration, therefore an equipment malfunction.

Environmental conditions during manufacture

15-30 °C (room temperature), RH 30-80 %.

Storage conditions

T = 10-30 °C, RH 30-80 %, keep out of direct light.

Transportation Conditions

Storage conditions should be strictly followed. Transport in closed and reinforced boxes in order to avoid damages. Transportation of this product does not represent a risk for health.

Shelf-life

Metal capsule: 500 cycles. Plastic helix: 250 cycles (replaceable at very low cost).

Packing

1 Helix-PCD + 1 Cotton bag

Packing information: product code and description, manufacturer information and data on label.

Labelling

Only on product's box: product code and batch number.

Possible target markets

Healthcare, Food, Pharmaceutical and Medical Industries.

Other important information (optional)

Read product's directions for use thoroughly before use.

Note: When necessary, limits for each product's parameters should be specified.

TERRAGENE

Description of identified and/or predictable risks which could represent a risk for the product.

Do not place the Helix-PCD device within the cotton bag when sterilizing with H₂O₂ and Formaldehyde; since it can interfere with these processes.

Do not store the product near sterilizing agents.

Do not expose this product to chemical agents.

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