
IT10

Indicador Integrador fisicoquímico

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IT10

Physicochemical Integrator Indicator

For Ethylene Oxide sterilization processes

Usage

Control of Ethylene Oxide (EO) sterilization processes.

Applicable Regulation

EN ISO 11140-1:2009

Classification

Class 1, according to risk.

Authorization

Designed under Quality Management Systems standards ISO 13485:2003/NS-EN ISO 13485:2012. ANMAT (Argentinean National Administration of Drugs, Food and Medical Technology) PM 1614-4.

Characteristics

Class 5 Integrator Indicator.

White polypropylene base: 17.2 mm x 5.0 mm x 84.4 mm

Auto adhesive label: 81.6 mm

Glass capillary tube: 74 mm long x 1.6 mm external diameter, 1.0 mm internal diameter. Contains indicator resin.

Polyethylene bag: 70 µm thick, 40 mm wide, 110 mm long. Thermally sealed on both ends.

Initial color: **orange**

Final color: **blue**

IT10 integrator indicator is specifically designed to react with all critical parameters of EO sterilization cycle (temperature, gas concentration, humidity and sterilization time). The integrating system permits to establish a gradient concentration of sterilizing agent which advances through the capillary tube reacting with the packed resin. IT10 integrator indicator has a reference pointer on the label (orange triangle), which indicates the minimum acceptance limit. This is the minimum limit of color change for which the required critical parameter values are accepted to meet an inactivation status of a spore population with D value specified between the limits accepted by the corresponding standards. The reference pointer corresponds to the kill time of a 10^6 *Bacillus atrophaeus* ATCC 9372 spore population (Bionova® biological indicator), calculated in BIER (Biological Indicator Evaluator Resistometer).

Conditions: EO 600 mg/l, RH 60 %, 54 °C // EO 600 mg/l, RH 60 %, 37 °C.

Environmental conditions during manufacture

15-30 °C (room temperature), RH 30-80 %.

Storage conditions

T = 10-30°C, RH 30-80 %, keep out of direct light.

Transportation conditions

Storage conditions should be strictly followed. Transport in closed and reinforced boxes in order to avoid damages. Transportation of this product does not represent a risk for health.

Shelf-life

5 years

Packing

100 units per box.

Packing information: product code and description, process for intended use, operating instructions, product presentation and manufacturer information.

Labelling

On product: product description, classification according to regulation, process for intended use, batch number and expiration date.

On product's box: product description, batch number and expiration date.

Possible target markets

Healthcare, Food, Pharmaceutical and Medical Industries.

Other important information (optional)

Read product's directions for use thoroughly before use.

Note: When necessary, limits for each product's parameters should be specified.

Description of identified and/or predictable risks which could represent a risk for the product

Do not store the product near sterilizing agents.

Do not expose this product to Steam, Dry Heat or any sterilization processes other than EO.

Do not reuse the sterilizer until the indicator has changed into the color indicated in the product.