# IT27

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
# IT27

## Three Level Emulator Indicator

*For Steam sterilization processes.*

### Usage

Control of Steam sterilization processes.

**Emulation conditions:** 7 minutes, 134 °C // 20 minutes, 121 °C.

### Applicable Regulation

ISO 11140-1:2014

### Classification

Class 1, according to risk.

### Authorization

Designed under Quality Management Systems standards ISO 13485:2003/NS-EN ISO 13485:2012. ANMAT (Argentinean National Administration of Drugs, Food and Medical Technology) PM 1614-4.

### Characteristics

Type 6 Emulator Indicator.

One-side-laminated card printed with indicator ink.

Size: 60 mm x 25 mm

**Level 1 – Blue:** Exposure level, it turns to black to indicate the exposure to the process but it does not assure that the conditions indicated in the emulator have been met (UNSAFE condition).

**Level 2 – Green:** Emulation level, it turns to black thus indicating that conditions indicated in the emulator have been met (SAFE condition). The SAFE condition is reached in steam sterilization processes of 20 minutes at 121 °C or 7 minutes at 134 °C.

**Level 3 – Yellow:** Over-emulation level, it turns to black when emulation conditions indicated in the product have been exceeded (EXTENDED condition).

### Environmental conditions during manufacture

15-30 °C (room temperature), RH 30-80 %.

### Storage conditions

T = 10-30 °C, RH 30-80 %, keep out of direct light.

### Transportation conditions

Storage conditions should be strictly followed. Transport in closed and reinforced boxes in order to avoid damages. Transportation of this product does not represent a risk for health.

### Shelf-life

5 years

**Packing**

250 units per aluminium foil bag.

Packing information: product code and description, process for intended use, operating instructions, product presentation, classification according to regulation, batch number, manufacturer information and expiration date.

**Labelling**

Information only on product's bag.

**Possible target markets**

Healthcare, Food, Pharmaceutical and Medical Industries.

**Other important information (optional)**

Read product's directions for use thoroughly before use.

**Note:** When necessary, limits for each product's parameters should be specified.

**Description of identified and/or predictable risks which could represent a risk for the product**

Do not store the product near sterilizing agents.

Do not expose this product to EO, Dry Heat or any sterilization processes other than Steam.

Do not reuse the sterilizer until the indicator has changed into the color indicated in the product.