# KCCDER2

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# KCCDER2

## Cleaning Challenge Device Kit

*For Automatic Endoscope Reprocessors and/or Automatic Thermidisinfecting Washers.*

### Usage

Cleaning Challenge Device Kit designed to evaluate the effectiveness of flexible endoscopes cleaning or other instruments to be reprocessed in Automatic Endoscope Reprocessors and/or Automatic Thermidisinfecting Washers.

### Applicable Regulation

ISO/TS 15883-5:2005

### Classification

Class 1, according to risk.

### Authorization

Designed under Quality Management System standard ISO 13485:2003/NS-EN ISO 13485:2012.

### Characteristics

Challenge Device

Chemdye® CCDER2 Kit consists of a perforated plastic capsule, which is able to hold CDWE cleaning indicator. The capsule is connected to a high resistance plastic tube of 2 m and 2 mm internal diameter. Also, the tube has a luer-lock connector, which enables connection to Automatic Endoscope Reprocessors (AERs) and/or Automatic Thermidisinfecting Washers with available connector terminals.

CDWE

Synthetic, self-adhesive and temperature-stable support printed with indicator colored ink. Allows detection of all factors affecting the outcome of the cleaning process.

Size: 38 x 6 mm; 20 strips per sheet.

Sheet size: 98.5 mm x 80 mm.

Initial Color: fuchsia

Result: elimination of fuchsia organic substance.

The indicator ink is non-toxic

NOTE: If unfavorable results are obtained, check quality/quantity of chemicals, water pressure, temperature, cleaning time and general condition of the washing machine.

### Environmental conditions during manufacture

15-30 °C (room temperature), RH 30-80%.

### Storage conditions

T = 10-30 °C, RH 30-80%, store in dark place.

### Transportation conditions

Storage conditions should be strictly followed.

Transport in closed and strengthened boxes in order to avoid damages. Transportation of this product does not represent a risk for health.

### Shelf-life

Plastic capsule: 250 cycles.

Cleaning indicator: 5 years

**Packing**

1 Challenge Device + 120 indicator strips + 1 Luer-lock connector  
Packing information: product code and description, intended use, classification according to regulation, operative instructions, manufacturer information and data on label.

**Labelling**

On product's box: Product code and description, intended use, classification according to regulation, batch number and expiration date.

**Possible target markets**

Healthcare and Medical Industries.

**Other important information (optional)**

Read product's directions for use thoroughly before use.

**Description of identified and/or predictable risks which could represent a risk for the product**

Store the product within its' original package, away from light, disinfectants, gases and humidity.

Do not re-use the indicator.

Do not touch indicator ink with your hands.