KH2X15-7.0Y/P

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KH2X15-7.0Y/P

Helix-PCD + Chemical Indicator test strips Kit.

For Steam sterilization processes control.

Usage

Process Challenge Device for Steam sterilization processes control.
For safe release hollow loads.
Conditions: 7 min at 134°C

Applicable regulation

EN 867-5, EN 285 and ISO 11140-1:2014

Authorization

Designed under Quality Management System standards ISO 9001:2008 and ISO 13485:2012.

Characteristics

Helix-PCD

Plastic capsule of great resistance and durability, to hold a Chemical Indicator designed for this system. The capsule is connected to a high resistance plastic tube of 1.5 m length (see catalogue for other lengths) and 2 mm internal diameter. Special silicone connector coated by external connector.

PCD-A-7.0Y

Type 2 Chemical Indicator.

Self-adhesive strip printed with indicator ink. Size: 77 x 6 mm.
10 strips per sheet. Size: 80 x 100 mm.

Initial color: yellow

Final color: black

NOTE: If the 4 bars change to the final color, the sterilization cycle was successful, the load can be securely released. If one or more bars did not turn to the final color, this shows a failure in sterilizing agent penetration, the load cannot be considered sterile, that is, it cannot be released.

Environmental conditions during manufacture

15-30 °C (room temperature), RH 30-80 %.

Storage conditions

T = 10-30 °C, RH 30-80 %, store in dark place.

Transportation conditions

Storage conditions should be strictly followed.
Carry in closed and reinforced boxes in order to avoid damages.
Transportation of this product does not represent a risk for health.

Shelf-life

Helix-PCD: 250 cycles, PCD-A-7.0Y: 5 years

Packing

100 or 250 PCD-A-7.0Ystrips in sheets of 10 units + 1 Helix-PCD + 1 cotton bag.

Packing information: product code and description, process for intended use, presentation, manufacturer information and data on label.

Labelling

On product's box: product code and batch number.

On helix-PCD bag: product code and batch number.

On PCD-A-7.0Ybox: batch number and expiration date.

Possible target markets

Healthcare, Food, Pharmaceutical and Medical industries.

Other important information (optional)

Read product's directions for use thoroughly before use.

Description of identified and/or predictable risks which could represent a risk for the product

Do not store the product near sterilizing agents.

Do not expose this product to EO, Dry Heat, Formaldehyde or any sterilization processes other than Steam.

Do not reuse the sterilizer until the color change of the chemical indicator is the one indicated in the product.