Auto-reader

for Rapid and Super Rapid
Biological Indicators
and Hygiene Monitoring
systems

Innovation in technologies for sterilization and disinfection control


IC10/20FR

TERRAGENE®

Rev.22 | June 2017

Constituent parts

Description of use

Bionova® IC10/20FR incubator has been designed to incubate BT110 Rapid Readout Self-Contained Biological Indicators (Bls) for Ethylene Oxide sterilization processes control, BT220 and BT221 for Steam sterilization processes, BT102 for Formaldehyde sterilization processes control, BT95 for Plasma or Vaporized Hydrogen Peroxide sterilization and Super Rapid Readout Self-Contained Biological Indicators BT222 and BT223 for Steam sterilization processes control.

Microorganisms contained within the Biological Indicators vary depending on the sterilization process for which they have been designed, therefore, incubation temperatures differ depending on the indicator. Indicators for Ethylene Oxide (EO) sterilization control, should be incubated at 37°C, while those involved in the monitoring of the rest of the sterilization processes should be incubated at 60°C.

! Simultaneous incubation of Bls for EO processes control within Bls monitoring the other processes, is not possible.

Bionova® IC10/20FR incubator allows easy and rapid detection of positive and negative Biological Indicators. Five incubation programs are available: 4 hours at 37°C, 3 hours at 60 °C, 2 hours at 60 °C, 1 hour at 60 °C and 30 minutes at 60 °C. A positive result can also be evidenced through color change of the culture medium, which turns to yellow after 48 hours or 7 days of incubation at 60°C for BT220/BT221/BT222/BT223/BT102/BT95 and after a 7-day incubation at 37 °C for BT110 Bls. Incubation times depend on internal protocols of the laboratory or hospital. The incubator allows selection of incubation time programs for two zones in an independent way. Zone 1 includes positions 1 to 6, and zone 2 comprises positions 7 to 12.

The IC10/20FR Bionova® incubator contains a special cavity to run quantitative readings of the Chemdye® PRO1 System, destined to surface cleaning and contamination control.

Bionova® IC10/20FR auto-reader not only allows the detec-

tion of results, but also provides a printed ticket to record them. In addition, the auto-reader can be connected to computers by USB. This allows to link the results with the Bionova® Automatic readout and traceability system for sterilization biological indicators program. To download the latest software version go to: <http://www.terragine.com.ar/products/infection-control/control-of-sterilization-processes/incubators/ic1020f>, you can also download its operator's manual from there.

- 1 Protective cover.
- 2 Incubation area for 12 Bls. Protein pen detection. Bls Ampoule crusher.
- 3 Control Panel.
- 4 Red light/Positive indicator/Position Light.
- 5 Green Light/Negative indicator/Position Light.
- 6 Reading position number.
- 7 Program selector.
- 8 Temperature Selector.
- 9 Alarm Cancellation Button.
- 10 Thermal printer.
- 11 Printer's paper indicator light.
- 12 Paper's traction button.
- 13 Cavity for paper.
- 14 USB port.
- 15 Input for power source plug (12V DC).
- 16 Hole for external temperature control.

Safety information

Symbols

Operating relative humidity.

Operating environmental temperature.

Caution, warning.

Caution: Risk of electric shock.

Caution: Hot surface.

Important, attention.

Manufacturer.

European conformity.

Serial number.

Direct current.

Product denomination.

To avoid risks and damaging the equipment

- Do not place the incubator in an environment exposed to direct sunlight or high intensity light lamps.
- Do not place the incubator near devices that emit electromagnetic fields.
- Do not use the equipment on vibrating surfaces.
- Do not pour any liquid inside.
- Do not immerse into any liquid.
- Use indoors only.
- Disconnect the power cord before cleaning.
- Do not use abrasive, corrosive cleaners or disinfectants.
- In case of technical fault, contact the manufacturer for support. Please, do not try to open or repair the incubator on your own, since this will imply the loss of product warranty and could lead to a major and irreversible damage.
- Make sure that the incubator is connected to a properly rated power cord.

To reduce the risk of using incompletely sterilized loads:

-Make sure to use each BI with its incubation and reading program:

BT110 (green top) with 4 hs. (37 °C) program

BT220 (brown top) with 3hs. (60 °C) program

BT221 (blue top) with 1 h. (60 °C) program

BT222 (violet top) with 1 h. (60 °C) program

BT223 (light blue top) with 30 min. (60 °C) program

BT102 (red top) with 2 hs. (60 °C) program

BT95 (white top) with 2 hs. (60 °C) program

-Check that culture medium completely wets the spore carrier.

-Do not remove the BI until the equipment reports the final readout result, check the impression of the ticket.

To avoid the risk of injury, because of glass fragments generated after the glass ampoule inside the BI tube is crushed:

-Cool the BI during the indicated time before crushing the ampoule.

-Do not handle the BI excessively since this might cause the glass ampoule to burst.

-Wear safety glasses and gloves when removing the BI from the sterilizer.

-Wear safety protective glasses and gloves to crush the ampoule and to press the cap of BI.

-Do not use your finger for crushing BI.

To avoid a potentially hazardous situation: avoid contact with the hot surface of the metal incubator block.

To avoid BIs from absorbing fluorescent residues:

avoid BI's direct contact with chemical indicators or tapes.

- Do not use this product in a manner not specified by Terragene® S.A., otherwise the protection provided by the product might be affected.

Do not replace the power supply cord for other not provided within the product.

Only personnel authorized by Terragene® S.A can access the inside of the incubator and its components. Parts or components inside the machine should not be manipulated by the user.

Safety information

Power supply specifications

Input parameters	Operating Condition	Units
Voltage range	(100-240)	Volts
Frequency	50/60	Hertz
Current	0,2	Ampers
Output parameters	Values	Units
Voltage	12	DC Volts
Current	4	Ampers

Terragene® recommends using a UPS instead of stabilizers, since the former has two functions: it stabilizes and provides a continual power source during energy cutoff.

Environment operating conditions

Environmental cond.	Functional cond.	Units.
Altitude	3500 (máx.)	Meters
Operation temperature	10-30	°C
Relative humidity	30-80	%
Installation/overvoltage	category II	
Degree of contamination	2	
Storage temperature	0 - 50	°C
Voltage	12	DC volts

Serial/Batch Number

For easy identification, each IC10/20FR incubator has a unique serial number and a batch number printed on a label found at the bottom of the incubator and on its packaging. Please, register your serial number and batch number in this manual for future reference;

Serial Number:

Batch Number: _____

Record Table

Use the following table to record firmware updates.

NOTE: manufacturer firmware is printed on a label found on the bottom of the incubator and on its packaging.

Please, register your manufactory firmware in this manual for future reference:

Manufacture Firmware: _____

Instructions for use

Start-up

-  All equipment is manufactured with the following settings:
- Time zone: UTC +0:00
- Printing language: English

For changing time zone, refer to *Time zone change mode* section.

For changing printing language, refer to *Printing language* section.

1-Place the auto-reader on a firm surface, free from vibrations, away from sunlight, currents of hot or cold air, chemicals and corrosive or flammable substances. Do not place the equipment in a way that disconnection of the plug from the power supply could be difficult. Leave a gap of at least 10 cm from the wall. Do not move the auto-reader periodically or during its use. Connect the auto-reader to a secure and stable electrical power source.

-  Do not wet or heat the auto-reader. If liquid is spilled on the incubator, disconnect it and dry it immediately. Before turning on, verify that all reading positions are empty.

2-Turn on the machine by connecting the end of the AC power supply to the mains and then connect the other end of the source (plug) on the rear of the IC10/20FR unit. The printer's paper indicator light will turn on, if this light becomes intermittent, this means that the thermal paper must be loaded, or that the printer cover is not tightly closed. Otherwise, the light will remain stationary.

The machine will start running the last time and temperature set-up.

-  If the incubator has an out-of-date Firmware version, the settings might not be saved.

3-Synchronize IC10/20FR with local time using Bionova® software available at <http://www.terragene.com.ar/productos/control-de-infecciones/control-de-procesos-de-esterilizacion/incubadoras/ic10-20fr/>.

For this, connect a PC to the USB plug of the incubator and start the program. The Bionova® IC10/20FR will be synchronized with the time in the PC. Synchronization can be done as long there is not an ongoing reading.

4- Check that the incubation temperature is correct. To change it, press the  button for 3 seconds. Lights corresponding to 37 °C and 60 °C will start blinking, indicating that a temperature program must be selected. For temperature selection, press the  button to select 37 °C or 60 °C. After 4 seconds the temperature modification will be accepted.

For 37 °C, the IC10/20FR auto-reader will automatically select the 4 hrs. final reading program in both halves of the incubation area (1-6 and 7-12).

For 60 °C, 30 min., 1, 2 or 3 hrs. program can be selected. It's possible to select different programs for each half or the same for both, depending on the indicator's final reading characteristics.

Check that the time setting is appropriate for the BIs to be incubated in each zone.

To change the time program, press the  button of the correspondent zone as many times as needed to choose the right setting.

5-Wait for temperature of the machine to be stable. Once is stable, the indicator of the desired temperature will remain steady.

-  Do not put any biological indicator in the reading positions until the incubator does not confirm temperature stability.

The auto-reader will keep the selected incubation program during the reading process. The system does not allow to modify the incubation program once the reading process has begun, unless the reading is finished.

6-Start readings once the blue temperature selection's light remains stable, this indicates that incubation temperature was reached. Once the selected temperature has been reached and remains stable, after turning on the equipment, a one-time automatic check of the Biological Indicators and Protein Pen incubation positions will be carried out. Indicating lights of each well will turn on bringing the result, a green light will show a satisfactory result while a red light will indicate error.

When the equipment detects an error in any position, a ticket indicating the day, time and position involved will be printed. To guarantee the reliability of the results, these positions will be deactivated, thus preventing their use. To evidence this, the red light of that position will remain blinking indefinitely, in order to indicate that the position has been disabled.

-  For a correct automatic check, the protective cover of the incubator should be in its place.

-  Read BIs instructions for use. Before placing the vials in the incubator, press the top to seal the tube. Crush the ampoule (using the ampoule crusher within the box of BIs or the one located within the incubation area of the incubator) and make sure that the liquid has completely wet the spore carrier placed at the bottom of the biological indicator plastic tube.

-  Do not remove or change placement of BI once the reading process has begun. If this occurs, results may be invalidated.

-  Use a non-sterilized biological indicator as a positive control every time a processed indicator is incubated. For more information refer to the *Positive Control* section of this user guide.

7-Place the rapid BI in a reading position and wait for the result. Reading will start automatically. Immediately after placing the BI in a position, an alarm will sound indicating that reading has successfully began and the position LED will blink until the result is informed.

Instructions for use

NOTE: If the incubator does not run automatically, it might indicate that the incubator temperature is not steady yet (check the temperature indicating light to remain steady) or that the incubation position is disabled due to an error in the initial autotest.

8-If the incubator detects a positive result, the red light and the audible alarm will turn on. This will indicate that the sterilization process to which the BI was subject, has failed. 30 seconds after removing the IB from the incubation position, the alarm and the red light will turn off automatically and the position will be available to start a new reading.

The negative BIs will be informed with a green light in the corresponding position. This will indicate that the sterilization process to which it was subject has been successful. The green light will remain on for up to 30 seconds after the BI is withdrawn from the reading position.

The incubator will inform results within 30 min. (60 °C), 1 h. (60 °C), 2 hs. (60 °C), 3 hs. (60 °C) or 4 hs. (37 °C), according to the program selected. These periods of incubation are the time-limits used by the equipment to inform the fluorescence final readout.

9-Every time a positive result is detected, it will be informed through the printing of the corresponding ticket and with an audible alarm. The alarm can be cancelled by pressing the button .

10-To register the tickets, press paper's traction button, take the paper and press down for the paper to be cut with the serrated edge of the printer.

NOTE: If the paper runs out, printer's blue light will start blinking (see *Replacement of the paper roll in the printer*). The incubator will save the last 12 results, for re-printing, follow the *Reprinting of the last results* instructions.

BIONOVA IC10/20FR RAPID READOUT INCUBATOR PROGRAM: XX/XX DATE: DD/MM/YY START TIME: HH:MM READOUT TIME: HH:MM SAMPLE TUBE: NN POSITIVE	BIONOVA IC10/20FR RAPID READOUT INCUBATOR PROGRAM: XX/XX DATE: DD/MM/YY START TIME: HH:MM READOUT TIME: HH:MM SAMPLE TUBE: NN CANCELED	BIONOVA IC10/20FR RAPID READOUT INCUBATOR PROGRAM: XX/XX DATE: DD/MM/YY START TIME: HH:MM READOUT TIME: HH:MM SAMPLE TUBE: NN NEGATIVE
POSITIVE	CANCELED	NEGATIVE

Thermal paper specifications

Recommended paper: JUJO AF50KSE3 or similar (order code ICTP).

Paper width	58 mm
Maximum paper thickness	60 g/m2
Maximum diameter size	50 mm

Replacement of the paper roll in the printer

 The printer uses thermal paper rolls.

1. To replace the paper roll, pull the handle of the printer's door.

2. Open the printer's cover and remove the spent roll.

3. Place the new paper roll with the outer side up.

4. Close the printer cover by pressing on the sides of the lid.

Reprinting of results


The incubator allows reprinting of the last 12 results by fluorescence readout and the last result of hygiene monitoring reading. For printing the results, press the .

Remaining incubation time

The BIONOVA IC10/20FR Auto-Reader also allows verifying the remaining incubation time in every reading position. Whenever a reading process is in course, press and hold the .

Instructions for use

ticket will then be printed containing the following information: reading position, incubation program and remaining incubation time.

NOTE: This feature is available starting from the IC10/20FR Rev. 2.5 Firmware version.

For the detailed procedure for updating the firmware, please refer to the *Firmware Update section* of this user guide.

Positive control

The positive control should belong to the same batch of the processed BI.

Use a non-sterilized biological indicator as a positive control every time a processed indicator is incubated.

Incubate the processed biological indicator and the indicator used as positive control as long as specified in the instructions.

Press the cap to seal the tube and crush the ampoule, make sure the media completely wets the carrier. Write a C (for *control*) and the date on the label. The positive control can be placed in any well of the Bionova® IC10/20FR incubator.

It is a good practice to use a positive control, this helps ensuring:

- Correct temperature is reached
- Viability of spores has not been altered due to improper storage temperature, humidity, or proximity to chemicals.
- Capability of the media to promote rapid growth.
- Proper functioning of Bionova® IC10/20FR incubator.

Disposal

Dispose BIs after use, according to your healthcare and safety regulations (see BIs instructions for use).

Interpretation of results

Rapid Readout

The indicator used as positive control should show a positive result (red light). The processed BI results are not valid until the biological indicator used as positive control shows a positive result (red light).

A processed BI with a positive result (red), indicates a failure in the sterilization process. A negative result, (green) after 30 min. (60 °C), 1 h. (60 °C), 2 hs. (60 °C), 3 hs. (60°C) 4 hs. (37 °C) according to the used program, indicates that the sterilization process was effective. Take action on any positive result immediately. A positive result indicates a sterilization process failure may have occurred. Determine the cause of the positive BI following current facility policies and procedures. Always retest the sterilizer and do not use the sterilizer for processing loads until 3 consecutive BI results are negative.

Visual Readout

The Rapid readout system allows the confirmation of results by color change after 48 hours or 7 days incubation analysis (if 7-day readout is made, an humidified incubator will be required to avoid media dry out). If the sterilization process was not successful, the growth media contained in the indicator will turn to yellow. If the sterilization process was successful, the media will remain unchanged. Color change to yellow of the growth media shows a failure in the sterilization process. If after the incubation time (see biological indicators directions for use), no color change is noticed in the processed indicators, then the result is negative (the sterilization process was successful).

The color of the growth media in the BI used as positive control must turn to yellow during the process of incubation for the results to be valid. Record the positive ones and discard them immediately as shown in the biological indicators instructions for use.

Instructions for use

Hygiene monitoring systems

Chemdye® PRO1MICRO reading system

The Chemdye® PRO1MICRO reading system consists of an incubation position that allows the quantitative analysis of the products destinated to evaluate cleaning and contamination control.

IC10/20FR can run a quantitative analysis of PRO1MICRO Hygiene monitoring system, adjusted by a BSA (Bovine Seric Albumin) calibration curve. This analysis allows carrying out the traceability of the surface cleaning control.

This method presents a Detection limit of 1 µg and a Sensibility of 0,3 µg. Nevertheless, for a better interpretation of results, these are informed with a resolution of 0,5 µg.

! Readings will be able to initiate once the incubation temperature has been selected and is stable. The lights corresponding to the selected temperature (37 °C or 60 °C) will blink until the temperature is reached, and then stay stable. PRO1MICRO readings must be performed at 60 °C.

Reading Process:

1-Press and hold the **Pro** button for 1 second to enter the reading program. A beep will sound. Then wait three seconds to confirm the selection. The incubator will emit two beeps indicating that a program has been selected.

2-Use the PRO1MICRO pen to take the sample from the surface to be analyzed. For this, follow directions described in Figure 1 (Also explained on PRO1MICRO's directions for use).

3-Once the solutions contained within the device have been correctly mixed, proceed with incubation (remember that the reading cone must not contain the swab). For this, insert the pen on the incubation/reading position (Figure 2). Make sure to place PRO1MICRO as it is indicated on the incubator label. When it is inserted in the correct position, it cannot be rotated (Figure 2).

4-Press and hold the **Pro** button for 3 seconds to initiate the incubation and reading process (15 minutes 60 °C). The lights **✓** and **✗** will start blinking.

! Do not move or remove the pen from its position during the reading process, this can cause erroneous readings or even cancellation of the same.


5-Once the incubation process is finished, an audible alarm informs the completion of the program. The samples are read and their protein value is calculated from a BSA (Bovine Seric Albumin) standard curve. A ticket is printed with the calculated protein concentration. If the protein concentration is below the limit of detection of the method, a green light **○** will appear. To inform the method has detected protein, a

red light **✗** will turn on. Also one of the five PRO panel lights will turn on to indicate an approximated concentration range calculated between 1 to 5 µg. If the calculated value is above 5 µg all lights of the panel will turn on.


Read the directions for use of each product and be aware of the sensibility limits.

Quantitative Results: For IC10/20FR-2.10 (March 2017) firmware version or higher (others versions could report different values).


Approximate amounts of BSA expressed in µg


Figure 1

Remove the swab from the device. Carefully apply 3 drops of moisturizer on the swab or surface to be tested.


Swab thoroughly the desired surface


Place swab back into the device.


Instructions for use


Place swab downwards firmly to activate.


Shake vigorously and collect solution in the readout cone.


Slide swab upwards.


Readout cone without swab.


Shake vigorously and collect solution in the readout cone.


Incubate and read.

Slide swab downwards.

Shake for 15 seconds.

Insert the pen on the incubation/reading position.

Slide swab upwards.

Instructions for use

Audible alarm

An audible alarm will sound every time a positive result is detected in the Bionova® IC10/20FR incubator. The alarm allows the user to immediately detect a positive result without the need to visually control the equipment. The alarm can be cancelled by pressing the button .

Cancelling a reading

When a BI is removed from its position during incubation process, an audible alarm will sound to indicate, the tube should be returned to the corresponding position.

If the indicator isn't re-introduced into the reading position after 10 seconds of being withdrawn, the reading will cancel and the cancellation will be informed through the printing of a ticket.

Temperature monitoring

The Bionova® IC10/20FR incubator contains an internal temperature control. If the temperature falls outside of the specified range (37 ± 2 or 60 ± 2 °C), the blue temperature led light will start to twinkle.

Temperature can be externally monitored by using a Bionova® TB-IC1020 thermometer or other similar. The Auto-Reader IC10/20FR has a special hole for placing the thermometer which is located at the back of the device.

Care and cleaning

Clean the exterior surface of the Bionova® IC10/20FR incubator with a damp cloth moistened with an small amount of detergent. Do not clean interior components. If additional cleaning of the wells is required, contact your distributor or manufacturer.

 Always unplug the IC10/20FR incubator and allow it to cool before cleaning. Do not immerse the unit in liquid.

NOTE: Bionova® IC10/20FR incubator does not need routine maintenance.

Firmware update


The IC10/20FR auto-reader allows updating the firmware periodically (program that controls the device and defines its different functionalities) through the use of the Bioupdate Software. This software verifies, by using the internet, the latest firmware version available for IC10/20FR auto-reader; downloads it and installs it in a few seconds without the loss of any device data. The Bioupdate Software must be downloaded, only once, from the following web site: <http://www.terragene.com.ar/productos/control-de-infecciones/con>

trol-de-procesos-de-esterilizacion/incubadoras/ic10-20fr/
Download file:  Software Bioupdate

UPDATE PROCEDURE:

Before starting the process, make sure that the Bionova® software is closed. Turn off the incubator for five seconds, then turn it on and follow this steps:

- 1- Connect the incubator to a PC using a USB cable and make sure to have an internet connection.
- 2- Open Bioupdate program.
- 3- Press Start button to update.

4- Wait until the software indicates that update process has been completed. The incubator will print a ticket indicating the new version of the Firmware.

NOTE: This process must be repeated every time that a new firmware version is available, this will be informed in the website.

 Do not modify custom parameters or use the "Full Flash" functionality. These functions were developed to be used only by qualified technicians.

Calibration


Incubator's calibration procedure is described in the software manual attached ("Automatic reading and traceability system for biological indicators for sterilization"). It can also be downloaded from <http://www.terragene.com.ar/products/infection-control/control-of-sterilization-processes/incubators/ic1020f>

Time zone change mode

This function allows you to advance or delay the internal incubator clock in order to adopt a different time zone.


Enter Time zone change mode

De-energize the incubator.

Instructions for use

Press and hold both  buttons simultaneously and energize the equipment. A beep will be emitted to indicate access to *Time zone change*.


NOTE: when entering *Time zone change* mode a green light will turn on in one of the incubation/reading position. This number corresponds to the current incubator time.

Time zone change

To advance time, press  button as many times as hours you need to advance. To delay time press  button as many times as hours you need to delay.


NOTE: 37 °C light indicates a.m. 60 °C light indicates p.m.

Complete time zone change

To save changes press and hold  button (Figure 4) for three seconds. A beep will sound indicating that process has been completed successfully. A ticket with the new set time will be printed.

Figure 4

Cancel time zone change


To cancel time zone change process without saving the changes, disconnect power source.

Time setting mode


Time setting functionality allows you establish a specific time regardless of incubator current time.

Enter *Time setting mode*

De-energize the incubator.

Press and hold both  buttons and  button simultaneously and energize the equipment. A beep will be emitted to indicate access to *Time setting mode*.

Instructions for use

Hours setting


Lights corresponding to 1 to 10 incubation/reading positions indicate the digit value which is set. Lights from 1 to 9 indicate digit value from 1 to 9. Light 10 indicates 0 value.

In time setting mode, light 1 from PRO sector will turn on to indicate ten hour modification and incubation light corresponding to the 10th position will turn on to indicate setting started.

NOTE: 24-hour format.


 To change hours and minutes digits, use  and  buttons. To increase digit value press  ; to decrease press  (see Figure 5).

Figure 5

- Hours setting: Ten hours (H:h:m)

The blue light number 1 corresponding to PRO sector should be on. If any other light from PRO sector is on, press  button as many times as necessary to turn on the blue light number 1. This indicates ten hour setting.

To increase ten hours value, press  button. To reduce ten hours value, press  button. (See Figure 5).

NOTE: the possible values for ten hour are 0, 1 and 2.

- Hours setting: Unit hours (hH:m)

Press  button to turn on the blue light number 2. This indicates unit hours setting. If any other light from PRO sector is on, press  button as many times as necessary to turn on the blue light number 2.

To increase unit hours value, press  button. To decrease unit hours value, press  button. (See Figure 5).

Minutes setting

- Minutes setting: Ten minutes (hh:M)

Press  button to turn on the blue light number 3. This indicates ten minutes setting. If any other light from PRO sector is on, press  button as many times as necessary to turn on the blue light number 3.

To increase ten minutes value, press  button. To decrease ten minutes value, press  button. (See Figure 5).

NOTE: possible values for ten minutes are 0, 1, 2, 3, 4 y 5.

- Minutes setting: Unit minutes (hh:m)

Press  button to turn on the blue light number 4. This indicates unit minutes setting. If any other light from Pro sector is on, press  button as many times as necessary to turn on the blue light number 4.

To increase unit minutes value, press button. To decrease unit minutes value, press button. (See Figure 5).

Finish hours setting

To save changes press and hold button (see Figure 4) for three seconds. A beep will sound, PRO green light will turn on and a ticket with the new set time will be printed to indicate the procedure has been completed successfully.

NOTE: If the modification has not been done correctly, the incubator will emit 3 beeps and the red light from PRO sector will turn on. Repeat the procedure.

Cancel hours setting

To cancel hours setting procedure without saving the changes, disconnect power source.

Printing language change mode


This function allows you to change the incubator printing language.

NOTE: each language has a code number assigned:


- 1 → English
- 2 → Spanish
- 3 → Portuguese
- 4 → Turkish

Enter to Printing language change mode


De-energize the incubator.

Press and hold button and energize the equipment. A sound will be emitted to indicate Access to Printing language change mode.

Use and buttons to select the position number corresponding to language code.

Finish printing language change

To save changes press and hold button for three seconds. A sound will be emitted indicating that process has been completed successfully. A ticket with the new language will be printed.

Troubleshooting chart

Fault	Possible cause	Action
The incubator does not start	Power source is not connected.	Check that the power source is connected to the power supply and that the plug is connected to the incubator.
The incubator gives an error in one position during "Autotest".	A BI is placed in the position at the moment of starting the incubator.	Check that every position is empty at the moment of starting the incubator.
The incubator gives an error in the "Autotest".	Soil particles might be obstructing the light path between the sensors.	Ensure cleanliness of the incubation position. Use of air is recommended. Do not use solid objects that could damage internal components. Restart the incubator.
The incubator does not run a BI reading. The red light of that position is blinking.	The autotest gives an error in that position, which becomes disabled.	Make sure that the position is empty when starting the incubator. Restart the incubator.
The incubator does not run a BI readout in any position.	The incubation temperature is not stable.	Wait for temperature to be stable.
The incubator does not run Chemdyne® Pro readings. The red light of that position is blinking.	The autotest gives an error in the PRO position, which becomes disabled.	Make sure that the PRO position is empty when starting the incubator. Restart the incubator.
When initiating the PRO incubation program the machine makes a sound but it does not initiate.	The incubation temperature is not stable.	Wait for temperature to be stable.
The incubator does not allow to change the incubation program.	Ongoing reading.	Wait for the readings to finish.
The printer does not print and the printer indicating light blinks.	The printer cover is not tightly locked.	Check that the cover is tightly closed.
	Printer without paper.	Place a new paper roll in the right direction.
The printer releases unprinted paper.	Paper roll is placed in the wrong direction.	Check the paper roll direction.
The Bionova® traceability software does not detect the incubator.	The incubator is turned off or the USB cable is disconnected.	Check that the incubator is turned on, the USB cable is connected and the COM port is correctly installed and accepted by the operative system.
The update of the incubator fails.	The incubator is connected to the Bionova® traceability software.	Close the program, restart the incubator and try again.

Warranty

Terragene® S.A. guarantees both the quality of the product material components and the quality of its manufacturing process. Should any material or manufacturing faults be detected within the warranty period (1 year from the time of purchase), the only obligation of Terragene® S.A. will be product repairing or substitution.

Limitation of liability

Terragene® S.A. shall not be liable for any loss or damage that may result from the unsuitable use of the equipment, negligence or user's full responsibility.

Intellectual property

Bionova® IC10/20FR auto-reader incubator may affect the intellectual property of other manufactures in some countries, invading certain aspects of patents related to this product. Terragene® S.A. shall not be liable for any damages or injuries that may result from the marketing of this product in those countries where Terragene® S.A. does not market it directly.

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Terragene S.A.
Güemes 2879 (S2002JBM)
Rosario. Santa Fe. Argentina.

Index

Constituent parts.....	3
Description of use	3
Safety information	4
Symbols	4
Power supply specifications	5
Environment operating conditions	5
Serial/Batch Number	5
Record Table	5
Instructions for use.....	6
Start-up	6
Thermal paper specifications.....	7
Replacement of the paper roll in the printer.....	7
Reprinting of results	7
Remaining incubation time	8
Positive control.....	8
Disposal.....	8
Interpretation of results	8
Hygiene monitoring systems.....	9
Audible alarm	11
Canceling a reading.....	11
Temperature monitoring	11
Care and cleaning	11
Firmware update	11
Calibration.....	11
Time zone change mode	11
Time setting mode.....	12
Printing language change mode.....	14
Troubleshooting chart	15
Warranty.....	16
Limitation of liability	16
Intellectual property	16
Technical assistance	16

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