Process Challenge Device (PCD)

Quality certification

Product code: _____
LOT _____

Bionova SCBI code: _____
Geobacillus stearothermophilus ATCC 7953

LOT _____

Población / Population _____ UFC / CFU
121 °C

Valor D / D - value _____ min.

Tiempo sobrevida / Survival time _____ min.
 Survival time = $(\log_{10} \text{labeled population} - 2) \times \text{labeled D-value}$

Tiempo de muerte / Kill time _____ min.
 Kill time = $(\log_{10} \text{labeled population} + 4) \times \text{labeled D-value}$

135 °C

Valor D / D - value _____ min.


Tiempo sobrevida / Survival time _____ min.
 Survival time = $(\log_{10} \text{labeled population} - 2) \times \text{labeled D-value}$

Tiempo de muerte / Kill time _____ min.
 Kill time = $(\log_{10} \text{labeled population} + 4) \times \text{labeled D-value}$

Valor Z / Z-value _____ °C

Record Card code: _____
LOT _____

Temperature SV (Time)
121 °C _____ min.
128 °C _____ min.
135 °C _____ min.

According to ISO 11140-1:2014 standard

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Process Challenge Device

For a rapid and easy monitoring of Steam sterilization processes

EN

Composition

Each pack consists of a stack of porous cards holding a self-contained biological indicator (SCB) tube that contains a population of *Gedacillus stearothermophilus* ATCC 7953 spores soaked on a carrier as well as growth indicator medium contained in a glass ampoule. The pack also contains a Record Card with a self-adhesive Type 5 Chemical integrator indicator and, on the outside, a Type 1 process indicator.

Product description

PCD is a device designed to simulate a load to be sterilized, and to pose a challenge to the sterilization process. These devices are used to evaluate the effective performance of the process. The steam Process Challenge Device has been designed for a rapid and easy monitoring of Steam processes at 132/135 °C x 4 minutes and processes at 121 °C x 30 minutes. It detects the inadequate air removal and Steam penetration in Steam sterilizers. This device allows release of routine loads, specially implants. Furthermore, it allows to perform the routine monitoring and periodic validation of the sterilizers (after repair, installation, relocation).

The device consists of a disposable pre-assembled package as outlined in ANSI/AAMI ST79 and contains a Self-contained Biological Indicator (SCB) and a Type 5 Chemical integrator indicator printed on The Record Card. The chemical integrator on the Record Card indicates that sterilizing conditions were met.

Warning!

Do not use the PCD for monitoring Ethylene Oxide, Dry Heat, Chemical Vapors, Radiation or any sterilization process other than Steam.

Do not reuse the SCB.

Storage

Best stored in original box under the next conditions: temperature 10-30 °C, RH 30-80%.

Do not freeze.

Do not store near sterilizing agents or other chemical products.

Directions for use

1- Place the pack inside a normally loaded steam autoclave. Place it in areas that you consider a priori most inaccessible for the sterilizing agent (high temperature steam). Generally, the center of the load and areas on the bottom shelf near the door and over the drain are problematic.

2- Run sterilization cycle:

- 132/135 °C ≥ 4 minutes
- 121 °C ≥ 30 minutes

3- After the sterilization process has finished, open the sterilizer door, wait for 5 minutes and remove the test pack.

4- Check that the external process indicator color on the outside of the test pack has changed. Open the test pack, wait 5 minutes and remove the SCB. Allow it to cool down to room temperature.

PRECAUTION! Wear safety glasses and gloves when removing the biological indicator from the sterilized test pack. **WARNING!** Do not crush or handle the biological indicator excessively, since this might cause the glass ampoule to burst.

5- Check the chemical integrator on the Record Card for correct exposure. Color change to the reference confirms that the inside of the pack has been exposed to sufficient sterilization conditions. Otherwise, check the sterilization process.

6- Check the chemical indicator on the label of the biological indicator. A color change to brown confirms that the biological indicator has been exposed to steam. **IMPORTANT:** This color change does not indicate that the process was sufficient to achieve sterility.

7- Identify the Bionova® BT220 SCB by writing the sterilizer number (in case of having more than one sterilizer), load number, and processing date on the label. Fill out the required information on the Record Card.

8- Crush the glass ampoule contained in the SCB. Then shake the tube down vigorously, with movements similar to those performed to lower the temperature in a mercury thermometer, until the medium reaches the base of the tube and soaks the spores carrier entirely. Incubate at 60±2 °C in a compatible incubator. **IMPORTANT:** Use a non-sterilized biological indicator as a positive control whenever a processed biological indicator is incubated. The positive control ensures that appropriate incubation conditions are met. Both the positive control indicator and the processed indicator should belong to the same batch. The indicator used as positive control must yield a positive result.

9- Incubate the processed biological indicator and the indicator used as positive control for 3 hours at 60 °C to get the final fluorescence result.

A 48-hour readout by visual color change is optional to confirm the 3-hour result. A positive fluorescence result (or growth medium color change after the 48-hour optional incubation) of the processed biological indicator means a sterilization process failure has occurred. If a negative fluorescence result is obtained (or growth medium remains the original color after the optional incubation period) the sterilization process was satisfactory.

10- Wait for the final results. Fill out the SCB results and adhere the self-adhesive whole Record Card or alternatively, only the area containing the integrator indicator. **WARNING!** Do not use the sterilizer until the biological indicator test results are negative.

11- Discard the test pack and the SCB immediately.

Disposal

Dispose used biological indicators according to your country's healthcare and safety regulations. The positive biological indicator can be autoclaved at 121 °C for at least 20 minutes, or at 132 °C for 15 minutes in a gravity displacement steam sterilizer, or at 134 °C for 10 minutes in a vacuum assisted steam sterilizer.

For the territory of Brazil: The use of this product is regulated by art. 99 of RDC 15/2012.

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