Process Challenge Device (PCD)

Quality certification

Product code: _____

LOT: _____

Bionova SCBI code: _____

Geobacillus stearothermophilus ATCC 7953

LOT: _____

 _____ _____

Población / Population _____ UFC / CFU

121 °C

Valor D / D - value _____ min.

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Process Challenge Device

For a rapid and easy monitoring of Steam sterilization processes

EN

Composition

Each pack consists of a stack of porous cards holding a self-contained biological indicator (SCBI) tube that contains a population of *Geobacillus stearothermophilus* ATCC 7953 spores soaked on a carrier as well as growth indicator medium contained in a glass ampoule. The pack also contains a Class 5 moving front integrator, a self-adhesive Record Card and, on the outside, a Class 1 process indicator.

Product description

PCD is a device designed to simulate a load to be sterilized, and to pose a challenge to the sterilization process. These devices are used to evaluate the effective performance of the process. The steam Process Challenge Device has been designed for a rapid and easy monitoring of Steam sterilization processes at 132/135 °C \geq 4 minutes and processes at 121 °C \geq 30 minutes. It detects the inadequate air removal and Steam penetration in Steam sterilizers. This device allows release of routine loads, specially implants. Furthermore, it allows to perform the routine monitoring and periodic validation of the sterilizers (after repair, installation, relocation).

The device consists of a disposable pre-assembled package as outlined in ANSI/AAMI ST79 and contains a Self-contained Biological Indicator (SCBI) and a Class 5 moving front integrator and a self-adhesive Record Card. The chemical integrator indicates that sterilizing conditions were met.

Warning!

Do not use the PCD for monitoring Ethylene Oxide, Dry Heat, Chemical Vapors, Radiation or any sterilization process other than Steam.

Do not reuse the SCBI.

Storage

Store in a dark place at temperatures between 10-30 °C, 35-70 % Relative Humidity.

Do not freeze.

Do not store near sterilizing agents or other chemical products.

Directions for use

1- Place the pack inside a normally loaded steam autoclave. Place it in areas that you consider a priori most inaccessible for the sterilizing agent (high temperature steam). Generally, the center of the load and areas on the bottom shelf near the door and over the drain are problematic.

2- Run the sterilization cycle:

- 132/135 °C \geq 4 minutes

- 121 °C \geq 30 minutes

3- After the sterilization process has finished, open the sterilizer door, wait for 5 minutes and remove the test pack.

4- Check that the external process indicator color on the outside of the test pack has changed. Open the test pack, wait 5 minutes and remove the SCBI. Allow it to cool down to room temperature.

PRECAUTION! Wear safety glasses and gloves when removing the biological indicator from the sterilized test pack. **WARNING!** Do not crush or handle the biological indicator excessively, since this might cause the glass ampoule to burst.

5- Check IT26-C moving front integrator for correct exposure. If the dark bar has reached the ACCEPT zone, this confirms that the inside of the pack has been exposed to sufficient sterilization conditions. Otherwise, check the sterilization process.

6- Check the chemical indicator on the label of the biological indicator. A color change to brown confirms that the biological indicator has been exposed to steam. **IMPORTANT:** This color change does not indicate that the process was sufficient to achieve sterility.

7- Identify the Bionova® BT220 SCBI by writing the sterilizer number (in case of having more than one), load number, and processing date on the label. Fill out the required information on the Record Card.

8- Crush the glass ampoule contained in the SCBI. Then shake the tube down vigorously, with movements similar to those performed to lower the temperature in a mercury thermometer, until the medium reaches the base of the tube and soaks the spores carrier entirely. Incubate at 60±2 °C in a compatible incubator. **IMPORTANT:** Use a non-sterilized biological indicator as a positive control whenever a processed biological indicator is incubated. The positive control ensures that appropriate incubation conditions are met. Both the positive control indicator and the processed indicator should belong to the same batch. The indicator used as positive control must yield a positive result.

9- Incubate the processed biological indicator and the indicator used as positive control for 3 hours at 60 °C to get the final fluorescence result.

A 48-hour readout by visual color change is optional to confirm the 3-hour result. A positive fluorescence result (or growth medium color change after the 48-hour optional incubation) of the processed biological indicator means a sterilization process failure has occurred. If a negative fluorescence result is obtained (or growth medium remains the original color after the optional incubation period) the sterilization process was satisfactory.

10- Wait for the final results. Fill out the SCBI and the IT26-C Integrator results and adhere the self-adhesive Record Card. **WARNING!** Do not use the sterilizer until the biological indicator test results are negative.

11- Discard the test pack and the SCBI immediately.

Disposal

Discard biological indicators after use according to your country's healthcare and safety regulations. The positive biological indicator can be autoclaved at 121 °C for at least 20 minutes, or at 132 °C for 5 minutes in a gravity displacement steam sterilizer, or at 134 °C for 10 minutes in a vacuum-assisted steam sterilizer.

For the territory of Brazil: The use of this product is regulated by art. 99 of RDC 15/2012.

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Fabricado por Terragene S.A. - Güemes 2879 - (2000)