# BT222-T

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# BT222-T

## Self-Contained Biological Indicator.

*Super Rapid Readout Fluorescence System.*

### Usage

Monitoring dynamic-air-removal (pre-vacuum) steam sterilization processes at 121 °C – 135 °C and gravity steam sterilization processes at 121 °C – 135 °C.

### Applicable Regulation

ISO 11138-1:2006 and ISO 11138-3:2006; IRAM 37102-1:1999 and IRAM 37102-3:1999

### Classification

Class 1, according to risk.

### Authorization

Designed under Quality Management System standards ISO 13485:2003/NS-EN ISO 13485:2012. ANMAT (Argentinean National Administration of Drugs, Food and Medical Technology) PM 1614-1.

### Characteristics

Polypropylene tube: 50.0 mm high x 8.5 mm external diameter. Wall thickness: 0.5 mm

Polypropylene microfiber cap at the bottom of the tube.

Polypropylene cap: 10.5 mm external diameter, 16.5 mm high. Wall thickness: 1.0 mm

Glass ampoule: 35.0 – 38.0 mm high. mm. External diameter: 6.5 mm. Culture medium 0.5 – 0.7 ml, purple color.

Cap filter: medical grade paper, 60 g/m<sup>2</sup>, 16.0 mm diameter.

Spore carrier: filter paper: 5.0 x 20.0 mm, 160 g/m<sup>2</sup>.

≥ 10<sup>6</sup> *Geobacillus stearothermophilus* ATCC 7953 spores per vial.

Each tube is individually identified by a serial number printed on its label.

Final fluorescence reading is performed after 1 hour-incubation at 60 °C (sensitivity: 97 %).

An optional visual pH color change confirmation could be made after 48 hours of incubation. If sterilization process has not been successful, culture medium will change to a greenish color first, and then to yellow during incubation at 60 °C, thus showing the presence of living spores. If sterilization process is successful culture medium will remain purple after the incubation process.

7-day readout is optional and not intended to be routinely performed; it is an initial validation of the 1 hour-reading. Fluorescence results may be compared to the 7-day visual reading.

NOTE: if 7-day readout is performed, a humidified environment will be required to avoid medium to dry out.

D-Value: Not lower than 1.5 minutes at 121 °C. Another D-value is estimated at 135 °C for the calculation of Z-Value.

Z-Value: Not lower than 6 °C.

### Environmental conditions during manufacture

Room temperature 15-30 °C, RH 30-80 %. Sterility conditions are necessary only during inoculation process performed in laminar flow.

### Storage conditions

T = 10-30 °C, RH 30-80 %, keep in a dark place.

### **Transportation conditions**

Storage conditions should be strictly followed. Products should be transported in closed and reinforced boxes in order to avoid damages. Product transportation does not represent any risk for human health.

### **Shelf-life**

24 months.

### **Packing**

50 units per box.

Packing information: product code and description, process for intended use, presentation, classification according to regulation, manufacturer information and data on box's label.

Weight per box: 143,7 g.

NOTE: manufacture date is calculated by subtracting 24 months to the expiration date.

### **Labelling**

On product: 17.0 mm x 33.0 mm polypropylene label. Printed in black and purple. 1.8 mm x 10.0 mm chemical indicator line, printed with Steam reactive ink (Color change to brown). Graph showing final fluorescence reading time. Product code, batch number, serial number, process for intended use and bacterial strain.

On product's box: product code, batch number, bacterial load, manufacture and expiration date.

### **Possible target markets**

Healthcare and Industry.

### **Other important information (optional)**

Incubation at 60 °C must be performed in Bionova® IC10/20FR Auto-Reader incubator.

Read product's directions for use thoroughly before use.

NOTE: When necessary, limits of each product's parameters should be specified.

### **Description of identified and/or predictable risks which could represent a risk for the product**

Do not store the product near sterilizing agents.

Do not expose this product to EO, Dry Heat, Radiation or any sterilization process other than Steam.

The positive biological indicator can be sterilized in gravity steam sterilizers at 121 °C for at least 20 minutes or at 132 °C for at least 15 minutes or at 134 °C for at least 10 minutes in vacuum-assisted steam sterilizers.