# KPCD200-1

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# KPCD200-1

## Steam Test Pack PCD Kit

*Process Challenge Device for Steam sterilization processes.*

### Usage

Detects inadequate air removal and steam penetration in dynamic-air-removal (pre-vacuum) and gravity displacement steam autoclaves.

**Condition:** 132/134 °C in dynamic-air-removal steam sterilization processes  $\geq 4$  minutes and 121 °C in gravity displacement steam sterilization processes  $\geq 30$  minutes.

### Applicable Regulation

ISO 11140-1:2014, ISO 11138-1:2006, ISO 11138-3:2006, ANSI/AAMI ST79

### Classification

Class 1, according to risk.

### Authorization

Designed under Quality Management System standards ISO 13485:2003/NS-EN ISO 13485:2012. ANMAT (Argentinean National Administration of Drugs, Food and Medical Technology) PM 1614-4.

### Characteristics

Porous cards system holding a Self-contained Biological Indicator and a Type 5 Integrator Indicator printed on PCDBI-1-RC self-adhesive Record Card, where data on sterilization cycle may be written. Two possibilities for documentation record: the whole Record Card can be adhered, or alternatively, only the area containing the indicator may be recorded.

Cardboard box (10.5 cm x 6.6 cm x 1.8 cm) with a self-adhesive Type 1 Process Indicator.

#### Process Indicator:

Initial color: **blue**

Final color: **dark grey or black**

#### PCDBI-1-RC Integrator Indicator (within the challenge device):

Initial color: **yellow**

Final color: **black**

Integration condition is calibrated against the kill time of a  $10^6$  *G. stearothermophilus* ATCC 7953 spores population, calculated in BIER (Biological Indicator Evaluator Resistometer). Conditions: saturated steam at 121 °C, 128 °C, 135 °C.

#### Self-contained Biological Indicator (within the challenge device):

$\geq 10^5$  *G. stearothermophilus* ATCC 7953 spores per vial in paper carrier. Glass ampoule with growth indicator medium. Label with chemical indicator line printed with Steam reactive ink (Color change: pink to brown)

**Final results:** incubation for 8 hours at 60 °C (sensitivity 97%)

NOTE: If sterilization process was not successful, indicator medium will turn yellow after incubation, thus indicating the presence of living *G. stearothermophilus* spores. If sterilization process was successful, indicator medium will remain purple after incubation.

**Environmental conditions during manufacture**

Room temperature 15-30 °C, RH 30-80 %.

**Storage conditions**

T = 10-30 °C, RH 30-80 %, keep in a dark place.

**Transportation conditions**

Storage conditions should be strictly followed. Carry in closed and reinforced boxes in order to avoid damages. Transportation of this product does not represent a risk for health.

**Shelf-life**

24 months.

**Packing**

25 PCD200-1 + 25 BT200 Biological Indicators per box.

Packing information: product description and manufacturer information.

**Labelling**

On product's box: product code and description, process for intended use, batch number and expiration date.

**Possible target markets**

Healthcare, Food, Pharmaceutical and Medical Industries.

**Other important information (optional)**

Read product's directions for use thoroughly before use.

**Note:** When necessary, limits for each product's parameters should be specified.

**Description of identified and/or predictable risks which could represent a risk for the product**

Do not store the product near sterilizing agents.

Do not expose this product to EO, Dry Heat, Radiation or any sterilization process other than Steam.

The positive biological indicator can be autoclaved at 121 °C for at least 20 minutes, at 132 °C for at least 15 minutes in gravity displacement steam sterilization, or at 134 °C for at least 10 minutes in vacuum steam sterilizer.