BT97 BioSurf

Indicador Biológico

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BT97 BioSurf

Self-Contained Biological Indicator

Super Rapid Readout Fluorescence System.

Usage

Monitoring of airborne and surface disinfection by Hydrogen Peroxide (H₂O₂).

Authorization

Designed under Quality Management System standards ISO 13485:2003/NS-EN ISO 13485:2012.

Characteristics

Tubes A and B: polypropylene, 45 mm high x 8.5 mm external diameter. Wall thickness: 0.5 mm

Special filter in tube B: polypropylene micro fibers of 17 mm diameter.

Cap A: polypropylene, 10.5 mm external diameter, 13 mm high. Wall thickness: 1.0 mm. This cap contains a 28.0 mm stem with a slot in which the carrier is located.

Cap B: polypropylene, 10.0 mm external diameter, 16.5 mm high. Wall thickness: 1.0 mm.

Glass ampoule: 38.0 - 40.0 mm height. External diameter: 6.5 mm. Wall 0.2 - 0.3 mm thick.

Culture medium 0.65 - 0.75 ml, purple color.

Cap filter: polyethylene fibers, 160 g/m², 17 mm in diameter.

Carrier in tube A: stainless steel 7.0 mm x 20.0 mm x 0.8 mm.

≥ 10⁶ *Geobacillus stearothermophilus* ATCC 7953 spores per vial.

Final fluorescence reading is performed after 1-hour incubation at 60 °C (sensitivity: 97 %).

An optional visual pH color change confirmation could be performed after 48 hours of incubation. If disinfection process has not been successful, culture medium will change to yellow during incubation at 60 °C, thus showing the presence of living spores. If disinfection process is successful culture medium will remain purple after the incubation process.

7-day readout is optional and not intended to be routinely performed; it is an initial validation of the 1-hour reading. Fluorescence results may be compared to the 7-day visual reading.

NOTE: If 7-day readout is performed, a humidified environment will be required to prevent medium from drying out.

D-Value: 2 mg/l H₂O₂, 50 °C.

Environmental conditions during manufacture

Room temperature between 15-30 °C, relative humidity 30-80 %. Sterility conditions are necessary only during inoculation and assemble of the biological indicator processes, performed in laminar flow.

Storage conditions

Best stored in original box under the next conditions: temperature between 10-30 °C, relative humidity between 30-80 %.

Transportation conditions

Storage conditions should be strictly followed. Products should be transported in closed and reinforced boxes in order to avoid damages. Product transportation does not represent any risk for human health.

Shelf-life

24 months.

Packing

50 units per box. 25 tubes A + 25 tubes B.

Packing information: product code and description, process for intended use, quantity, manufacturer information and data on box's label.

Packing information: product description, process for intended use, batch code, date of manufactured and expiry.

Note: manufacture date is calculated by subtracting 24 months to the expiration date.

Labelling

On product: 17.0 mm x 33.0 mm polypropylene label on each tube. On tube A there is a graph showing the final readout time by fluorescence, product code, batch number, process for intended use and bacterial strain, printed in black. The tube B contains the same information and a 1.5 mm chemical indicator line, printed with Hydrogen Peroxide reactive ink, which turns to green after exposure.

On product's box: product code, batch number, bacterial load, manufacture and expiration date.

Possible target markets

Healthcare and Industry.

Other important information (optional)

It is advisable to incubate at 60 °C in Bionova® IC10/20FR, Bionova® IC10/20FR LCD or Bionova® Mini Bio Auto-Readers incubators.

Read product's directions for use thoroughly before use.

Note: When necessary, limits of each product's parameters should be specified.

Description of identified and/or predictable risks which could represent a risk for the product

Do not store the product near sterilizing agents.

Do not expose this product to EO, Dry heat, Radiation sterilization processes or any disinfection process other than Hydrogen Peroxide.

The positive biological indicator can be autoclaved at 121 °C for at least 20 minutes, at 132 °C for at least 15 minutes in gravity displacement steam sterilization, or at 134 °C for at least 10 minutes in vacuum steam sterilizer.