# Touchscreen Auto-reader for Rapid, Super Rapid, Ultra Rapid Biological Indicators and Hygiene Monitoring systems

*Innovation in technologies for sterilization and disinfection control*


IC10/20FRLCD

TERRAGENE®

Rev.2 | December 2017

# Composition

## References

- 1 Protective Cover.
- 2 Incubation area for 12 Bls. Protein Pen Detection.  
Bls ampoule crusher.
- 3 Touch Screen Display
- 4 Thermal Printer
- 5 Indicator Light for Paper Replacement  
and Printer Functioning.
- 6 Paper's Traction Button.
- 7 Cavity for Paper.
- 8 12 volt DC Input.
- 9 Hole for External Temperature Control.
- 10 Ethernet port (RJ45).

# Composition

## Main screen

## Screen References

	PRO icon
	Positioning icons
	Status position light
	Date and Time 13:34 hs. 25/02/16
	Remaining time icon
	Alarm icon
	Temperature icon
	Start PRO icon
	Cancel PRO icon
	Program configuration icon
	Settings for all positions icon
	Results history icon
	Temperature stability lights 37 60
	Bionova® icon
	Network connection indicator / Bionova connection indicator

## Description of use

Bionova® IC10/20FRLCD Auto-reader has been designed for the incubation of the following Self-Contained Biological Indicators (IBs):

- Rapid Readout Self-Contained Biological Indicators: BT110 for Ethylene Oxide sterilization control, BT220 and BT221 for Steam sterilization control, BT102 for Formaldehyde sterilization control and BT95 for Plasma or Vaporized Hydrogen Peroxide sterilization control.
- Super Rapid Readout Self-Contained Biological Indicators: BT222 and BT223 for Steam sterilization control and BT96 for the control of Plasma or Vaporized Hydrogen Peroxide sterilization control.

- Ultra Rapid Readout Self-Contained Biological Indicators: BT224 for Steam sterilization control.

Microorganisms contained within the Biological Indicators vary depending on the sterilization process for which they have been designed, therefore, incubation temperatures differ depending on the indicator. Indicators for Ethylene Oxide (EO) sterilization control, should be incubated at 37 °C, while those involved in the monitoring of the rest of the sterilization processes should be incubated at 60 °C.

Simultaneous incubation of IBs for EO processes control within IBs monitoring the other processes, is not possible.

Bionova® IC10/20FRLCD Auto-reader allows easy and rapid detection of positive and negative Biological Indicators. Six incubation programs are available: 4 hours at 37 °C, 3 hours at 60 °C, 2 hours at 60 °C, 1 hour at 60 °C, 30 minutes at 60 °C and 20 min. at 60 °C.

A positive result can also be evidenced by culture medium color change to yellow, read the instructions for use of each indicator for more information. The choice to perform or not the extended incubation to confirm by color change depends on the internal protocols of each laboratory or hospital.

The IC10/20FRLCD Bionova® Auto-reader contains a special cavity to run quantitative readings of the Chemdye® PRO System.

The equipment can be connected to a PC through an Ethernet port. This allows to link the results with the Bionova® Traceability Software for Automatic Readout and Traceability of Rapid, Super Rapid, Ultra Rapid Biological Indicators and Hygiene Monitoring System. To download the latest version of the software, go to: <http://www.terragene.com.ar/productos/infection-control/incubators/software-downloads/>.

# Safety information

## Symbols

 80 %  
30 % Operating relative humidity.

 30 °C  
10 °C Operating environmental temperature.

 Caution, warning.

 Caution: Risk of electric shock.

 Caution: Hot surface.

 Important, attention.

 Manufacturer.

 European conformity.

 Serial number.

 Direct current.

 Product denomination.

### To avoid risks and damaging the equipment

- Do not place the auto-reader in an environment exposed to direct sunlight or high intensity light lamps.
- Do not place the auto-reader near devices that emit electromagnetic fields.
- Do not use the equipment on vibrating surfaces.
- Do not pour any liquid inside.
- Do not immerse into any liquid.
- Use indoors only.
- Disconnect the power cord before cleaning.
- Do not use abrasive, corrosive cleaners or disinfectants.
- In case of technical fault, contact the manufacturer for support. Please, do not try to open or repair the auto-reader on your own, since this will imply the loss of product warranty and could lead to a major and irreversible damage.
- Make sure that the auto-reader is connected to a properly rated power cord.

### To reduce the risk of using incompletely sterilized loads:

- Make sure to use each BI with its incubation and reading program:
  - BT110** (green top) with 4 hs. (37 °C) program
  - BT220** (brown top) with 3hs. (60 °C) program
  - BT221** (blue top) with 1 h. (60 °C) program
  - BT222** (violet top) with 1 h. (60 °C) program
  - BT223** (light blue top) with 30 min. (60 °C) program
  - BT224** (aquamarine top) with 20 min. (60 °C) program
- BT102** (red top) with 2 hs. (60 °C) program
- BT95** (white top) with 2 hs. (60 °C) program
- BT96** (grey top) with 30 min. (60 °C) program.

-Check that culture medium completely wets the spore carrier.

-Do not remove the BI until the equipment reports the final readout result, check the impression of the ticket.

### To avoid the risk of injury, because of glass fragments generated after the glass ampoule inside the BI tube is crushed:

- Cool the BI during the indicated time before crushing the ampoule.
- Do not handle the BI excessively since this might cause the glass ampoule to burst.
- Wear safety glasses and gloves when removing the BI from the sterilizer.
- Wear safety protective glasses and gloves to crush the ampoule and to press the cap of BI.
- Do not use your finger for crushing BI.

### To avoid a potentially hazardous situation:

avoid contact with the hot surface of the metal auto-reader block.

### To avoid BIs from absorbing fluorescent residues:

avoid BIs direct contact with chemical indicators or tapes.

 Do not use this product in a manner not specified by Terragene® S.A., otherwise the protection provided by the product might be affected.

Do not replace the power supply cord for other not provided within the product.

Only personnel authorized by Terragene® S.A can access the inside of the auto-reader and its components. Parts or components inside the machine should not be manipulated by the user.

## Safety information

## Power supply specifications

<b>Input parameters</b>	<b>Operating Condition</b>	<b>Units</b>
Voltage range	(100-240)	Volts
Frequency	50/60	Hertz
Current	0,2	Ampers
<b>Output parameters</b>	<b>Values</b>	<b>Units</b>
Voltage	12	DC Volts
Current	4	Ampers

Terragene® recommends using a UPS instead of stabilizers, since the former has two functions: it stabilizes and provides a continual power source during energy cutoff.

## **Environment operating conditions**

<b>Environmental cond.</b>	<b>Functional cond.</b>	<b>Units.</b>
Altitude	3500 (máx.)	Meters
Operation temperature	10-30	°C
Relative humidity	30-80	%
Installation/overvoltage	category II	
Degree of contamination	2	
Storage temperature	0 - 50	°C
Voltage	12	DC volts

**Serial/Batch Number**

For easy identification, each IC10/20FRLCD Auto-reader has a unique serial number and a batch number printed on a label found at the bottom of the auto-reader and on its packaging. Please, register your serial number and batch number in this manual for future reference:

**Serial Number:**

**Batch Number:** \_\_\_\_\_

## Record Table

Use the following table to record firmware updates.

**NOTE:** manufactory firmware is printed on a label found on the bottom of the auto-reader and on its packaging.

Please, register your manufactory firmware in this manual for future reference:

## **Manufactory Firmware:**

# Instructions for use

## Start-up

- !** All equipment is manufactured with the following settings:
- Time zone: UTC +0:00
  - Printing language: English
- If you wish to change language or time zone, go to *Settings* located into *Website of the equipment*.

1-Place the auto-reader on a firm surface, free from vibrations, away from sunlight, currents of hot or cold air, chemicals and corrosive or flammable substances. Do not place the equipment in a way that disconnection of the plug from the power supply could be difficult. Leave a gap of at least 10 cm from the wall. Do not move the auto-reader periodically or during its use. Connect the auto-reader to a secure and stable electrical power source.

- !** Do not wet or heat the auto-reader. If liquid is spilled on the auto-reader, disconnect it and dry it immediately. Before turning on, verify that all reading positions are empty.

2-Turn on the machine by connecting the end of the AC power supply to the mains and then connect the other end of the source (plug) on the rear of de IC10/20FRLCD unit. The printer's paper light will turn on, if this light becomes intermittent, this means that thermal paper must be loaded or that the printer cover is not tightly closed. Otherwise, the light remains stationary.

3-Enter *Website of the equipment* and set language, time zone and clock (see *Settings* located into *Website of the equipment* section).

4-Select incubation program and desired temperature for each readout position.

- !** Bionova® IC10/20FRLCD allows the selection of different incubation programs for each one of the readout positions.

- !** It also allows to configure the same program for the 12 positions simultaneously. To do this, press the icon . It will change its appearance to when simultaneous configuration is possible. For this, no position must have a readout in progress.

The device will start by taking the last incubation program and temperature settings used. To change the incubation temperature, press icon for 3 seconds to select 37 °C or 60 °C as appropriate. A screen will be displayed informing that temperature has changed.

- !** The incubation temperature can be modified only when there is no reading program in progress.

To select an incubation program for each of the positions you must press the icon of the corresponding position, ex. . A screen will be displayed informing the current program of that position:

To change it, press the icon . A screen will be displayed with the incubation program options:

Select the desired option. This configuration will be stored for that temperature until the user resets it.

- !** A change of incubation program can be made only in the empty positions.

- !** The auto-reader will keep the selected temperature (37 °C or 60 °C) during the reading process. The system does not allow to modify the temperature once the reading process has begun, unless the reading is finished.

- !** Since there is only one incubation program for 37 °C temperature, it will be automatically selected for all positions when this temperature is selected. Once selected the temperature of 60 °C again, the auto-reader will remember the last incubation program used in each position for that temperature.

5-Once the selected temperature has been reached and remains stable, after turning on the equipment, a one-time automatic check of the Biological Indicators and Protein Pen incubation positions will be carried out. Indicating lights of each position will turn on bringing the result, a green light will show a satisfactory result while a red light will indicate error.

When the equipment detects an error in any position, a ticket indicating the day, time and position involved will be printed. These positions will be deactivated, thus preventing their use.

- !** For a correct automatic check, the protective cover of the auto-reader should be in its place and all the positions empty.

Read BIs instructions for use. Before placing the vials in the auto-reader, press the top to seal the tube. Crush the ampoule (using the ampoule crusher within the box of BIs or the one located within the incubation area of the auto-reader) and make sure that the liquid has completely wet the spore carrier placed at the bottom of the biological indicator plastic tube.

# Instructions for use

**!** Make sure that the reading position you want to use is configured with the incubation program corresponding to the biological indicator you want to analyze.

**!** Do not remove or change placement of BI once the reading process has begun. If this occurs, results may be invalidated.

**!** Use a non-sterilized biological indicator as a positive control every time a processed indicator is incubated. For more information refer to the *Positive Control* section of this user guide.

6-Place the rapid BI in a reading position and wait for the result. Reading will start automatically. Immediately after placing the BI in a position, the auto-reader IC10/20FRLCD will emit an audible alarm to indicate that the readout started and will show the information screen of the initiated position. Status position light will turn violet and will blink until the result is informed.

**NOTE:** If the auto-reader does not run automatically, it might indicate that the auto-reader temperature is not steady yet (check the temperature indicating light to remain steady) or that the incubation position is disabled due to an error in the initial autotest.

**!** When a BI is removed from its position during incubation process, an audible alarm will sound to indicate the tube should be returned to the corresponding position.

If the indicator isn't re-introduced into the reading position after 10 seconds of being withdrawn, the reading will be cancelled and the cancellation will be informed through the printing of a ticket.

7-When a positive result is detected, the red light and the audible alarm will turn on, this will indicate that the sterilization process to which the BI was subject, has failed. 30 seconds after removing the IB from the incubation position, the alarm and the red light will turn off automatically and the position will be available to start a new reading.

The negative BIs will be informed with a green light in the corresponding position, this will indicate that the sterilization process to which it was subject has been successful. The green light will remain on for up to 30 seconds after the BI is withdrawn from the reading position. The auto-reader will inform results within 20 min. (60 °C), 30 min. (60 °C), 1 h. (60 °C), 2 hs. (60 °C), 3 hs. (60 °C) or 4 hs. (37 °C), according to the program selected. These periods of incubation are the time-limits used by the equipment to inform the fluorescence final readout. Every time a positive result is detected, it will be informed through the printing of the corresponding ticket and with an audible alarm. The alarm can be cancelled by pressing the button .

8-To register the tickets, press the button for paper traction, take the paper and pull it up to cut it with the serrated edge of the printer.

**NOTE:** If the paper runs out, the light of paper repla-

cement will start blinking (see *Replacement of the paper roll in the printer*).

 BIONOVA  
RAPID READOUT INCUBATOR  
BIONOVA - IC10/20FRLCD  
INDICATOR: XXXX  
PROGRAM: XXXX  
DATE: XXXX  
START TIME: XXXX  
READOUT TIME: XXXX  
SAMPLE TUBE: XXXX  
POSITIVE

 BIONOVA  
RAPID READOUT INCUBATOR  
BIONOVA - IC10/20FRLCD  
INDICATOR: XXXX  
PROGRAM: XXXX  
DATE: XXXX  
START TIME: XXXX  
READOUT TIME: XXXX  
SAMPLE TUBE: XXXX  
CANCELED

 BIONOVA  
RAPID READOUT INCUBATOR  
BIONOVA - IC10/20FRLCD  
INDICATOR: XXXX  
PROGRAM: XXXX  
DATE: XXXX  
START TIME: XXXX  
READOUT TIME: XXXX  
SAMPLE TUBE: XXXX  
NEGATIVE

## Thermal paper specifications

Recommended paper: JUJO AF50KSE3 or similar (order code ICTP).

Paper width : 58 mm

Maximum paper thickness: 60 g/m<sup>2</sup>

Maximum diameter size: 50 mm

## Replacement of the paper roll

**!** The printer uses thermal paper rolls.

1. Pull the handle of the printer's door.

2. Open the printer's cover and remove the spent roll.

3. Place the new paper roll with the outer side up.

4. Close the printer cover by pressing on the sides of the lid.

# Instructions for use

## Remaining incubation time | Configured programs

Bionova® IC10/20FRLCD Auto-reader allows to verify remaining time of incubation and the configured programs of each reading position in a very simple way. To do this, press on the icon . Then, a new screen will open indicating the remaining readout time of each position that is performing an incubation. In positions where there is no readout initiated, the name of the selected program will be displayed. Press to return to the main screen or just wait one minute to return automatically.

## Positive control

**!** The positive control should belong to the same batch of the processed BI.

Use a non-sterilized biological indicator as a positive control every time a processed indicator is incubated.

Incubate the processed biological indicator and the indicator used as positive control as long as specified in the instructions.

Press the cap to seal the tube and crush the ampoule, make sure the media completely wets the carrier. Write a C (for control) and the date on the label. The positive control can be placed in any well of the Bionova® IC10/20FRLCD Auto-reader.

It is a good practice to use a positive control, this helps ensuring:

- Correct temperature is reached
- Viability of spores has not been altered due to improper storage temperature, humidity, or proximity to chemicals.
- Capability of the media to promote rapid growth and the generation of fluorescence.
- Proper functioning of Bionova® IC10/20FRLCD Auto-reader.

## Disposal

Dispose BIs after use, according to your healthcare and safety regulations (see BIs instructions for use).

## Interpretation of results

### Rapid Readout

The indicator used as positive control should show a positive result (red light ). The processed BI results are not valid until the biological indicator used as positive control shows a positive result .

A processed BI with a positive result indicates a failure in the sterilization process. A negative result (green light ) after 20 min. (60 °C), 30 min. (60 °C), 1 h. (60 °C), 2 hs. (60 °C), 3 hs. (60°C) 4 hs. (37 °C) according to the used program, indicates that the sterilization process was effective.

The presence of positive results indicates a sterilization process failure may have occurred and implies that the load must be reprocessed immediately. Facing these results, verify good practices for load preparation and usage of biological indicators. Check the volume of the load to be sterilized as well as the BI's position inside the load. In case the BI has been correctly used, it should be interpreted that it is likely that the sterilizer has failed and maintenance is recommended.

## Visual Readout

The Rapid readout system allows the confirmation of results by color change after 48 hours or 7 days incubation analysis (if 7-day readout is made, an humidified environment will be required to avoid media dry out). If the sterilization process was not successful, the culture medium in the indicator will turn to yellow, as shown in the BIs Instructions for Use. If the sterilization process was successful, the indicator medium will remain unchanged. Color change to yellow of the growth media shows a failure in the sterilization process. If after the incubation time (see biological indicators directions for use), no color change is noticed in the processed indicators, then the result is negative (the sterilization process was successful).

**!** The color of the growth media in the BI used as positive control must turn to yellow during the process of incubation for the results to be valid. Record the positive ones and discard them immediately as shown in the biological indicators instructions for use.

## Hygiene monitoring system

### Chemdye® PRO1 MICRO reading system

The Chemdye® PRO1 MICRO reading system consists of an incubation position that allows the quantitative analysis of the products destined to cleaning and contamination control.

IC10/20FRLCD device can perform a quantitative analysis of the PRO1 MICRO Hygiene Monitoring System, adjusted by a BSA (Bovine Seric Albumin) calibration curve. This analysis allows carrying out the traceability of the surface cleaning control.

This method presents a Detection limit of 1 µg and a Sensibility of 0,3 µg. Nevertheless, for a better interpretation of results, these are informed with a resolution of 0,5 µg.

**!** Readings will be able to initiate once the incubation temperature has been selected and is stable. The lights corresponding to the selected temperature (37 °C or 60 °C) will blink until the temperature is reached, and then stay stable. PRO1 MICRO readings must be performed at 60 °C.

## Reading process

1-Press icon . Then, the following screen will be displayed:

# Instructions for use

2-Use PRO1 MICRO pen to collect the sample of the surface to be analyzed. To do this, follow the instructions described in Figure 1 (Also explained in the instructions for use of PRO1 MICRO).

3-Once the solutions contained within the device have been correctly mixed, proceed with incubation (remember that the reading cone must not contain the swab). For this, insert the pen on the incubation/reading position (Figure 2). Make sure to place PRO-1MICRO as it is indicated on the auto-reader label. Insert the pencil making sure that it stops at the bottom of the reading position, if it is in the correct position, check that it cannot be rotated (Figure 2).


4-Press icon to start the incubation and reading process (15 minutes at 60 °C). The equipment will return to the main screen and the PRO status position light will remain blinking until the reading ends.

**!** If for any reason the reading needs to be cancelled before it ends, PRO position should be re-entered by pressing icon . Then, press icon for 3 seconds. The reading will be cancelled and its corresponding ticket will be printed. The equipment will display the main screen automatically.


5-Once the incubation process is over, an audible alarm indicates the completion of the program. The samples are read and the protein concentration is calculated by a calibration curve with BSA (Bovine Serum Albumin). Then, a ticket with the calculated protein concentration is printed. If the concentration is below the detection limit, a green light will turn on. On the opposite, a red light will turn on if the method detects the presence of protein.

To visualize the displayed protein value on the screen, press the icon .


Finally press to return to the main screen.


Swab thoroughly the desired surface.


Place swab back into the device.


Press the swab downwards firmly to activate.


Slide swab upwards.

Shake downwards vigorously until the solution turns to green.


Slide swab downwards.


Shake for 15 seconds.

## Figure 1


Remove swab from the device. Carefully apply 2 drops of moisturizer on the swab or surface to be tested.


# Instructions for use


Slide swab upwards.

Shake vigorously and collect solution in the readout cone.


Readout cone without swab

Incubate and read.

**Figure 2**

Insert the pen in the incubation/reading position.

## Results history visualization

IC10/20FRLCD Auto-reader automatically stores the last 208 reading results. You can access these results by pressing the icon . After this, a screen will be displayed to show the results ordered chronologically. The ticket of the desired result will be reprinted by pressing on its corresponding line.

You can also access the visualization of the results by accessing the website of the equipment (see *Website of the equipment* section).

## Audible alarm

An audible alarm will sound every time a positive result is detected by the Bionova® IC10/20FRLCD Auto-reader. The alarm allows the user to immediately detect a positive result without doing a visual control of the equipment. The alarm can be cancelled by pressing the icon .

## Temperature monitoring

The Bionova® IC10/20FRLCD Auto-reader contains an internal temperature control. If the temperature falls outside of the specified range (37±2 or 60±2 °C), the blue temperature led light will start to twinkle. Temperature can be externally monitored by using a Bionova® TB- IC1020 thermometer or other similar. IC10/20FRLCD Auto-reader has a special hole for placing the thermometer which is located at the back of the device.

## Care and cleaning

Clean the exterior surface of the Bionova® IC10/20FRLCD Auto-reader with a damp cloth moistened with a small amount of detergent. Do not clean interior components. If additional cleaning of the wells is required, contact your distributor or manufacturer.

Always unplug the IC10/20FRLCD Auto-reader and allow it to cool before cleaning. Do not immerse the unit in liquid.

**NOTE:** Bionova® IC10/20FRLCD does not need routine maintenance.

## Network cable specifications

The use of a category 5 UTP cable with RJ45 connectors for Ethernet networks is recommended.

## Firmware update

IC10/20FRLCD Auto-reader allows updating the firmware periodically (program that controls the device and defines its different functionalities) through the use of the Bioupdate Software. This software verifies, by using the internet, the latest firmware version available for IC10/20FRLCD Auto-reader; downloads it and installs it in a few seconds without the loss of any device data.

The Bioupdate Software must be downloaded, only once, from the following web site: <http://www.terragenecom.ar/productos/infection-control/incubators/software-downloads/>.

Download file: [Download Bioupdate program](#)


To ensure the proper functioning of the Bioupdate software, before starting the update procedure, make sure:

# Website of the equipment

- That the Bionova® software is closed. Turn off the computer for 5 seconds and turn it on again.
- Connect your PC and the auto-reader to the same network, preferably through the specified network cable. Avoid connections of the PC and the auto-reader to the network through WiFi.
- That the firewall and the PC Anti-Virus do not block the execution of the Bioupdate software.

## UPDATING PROCEDURE:

- 1-Connect the auto-reader to the same network that your PC is connected through a specified Ethernet cable and make sure you have an internet connection.
- 2-Verify that the connection indicator on the main screen changes from  to . This indicates that the auto-reader has a connection to the network. Enter the auto-reader Information screen by pressing .
- 2-Open the Bioupdate program.

- 3-Select the equipment you want to update, check that the IP is correct. Press the Start button to update.
- 4-Wait until the software indicates that the update process has been completed. The auto-reader will issue a ticket indicating the new Firmware version.

**NOTE:** This process must be repeated every time a new version of Firmware is available, which will be informed on the website.

**NOTE:** When it is necessary to update the graphical environment, the update process may take 30 minutes or more. This will be notified with a message to the user from the Bioupdate software.

## Calibration

The calibration procedure of the equipment is described in the attached software manual ("Bionova® Traceability Software for Automatic Reading and Traceability of Rapid, Super Rapid and Ultra Rapid Biological Indicators and Hygiene Monitoring System"). You can also download it by going to <http://www.terragine.com.ar/productos/infection-control/incubators/software-downloads/>. You can check that the equipment is connected to the Traceability software if the Connection Indicator changes from  to .

The IC10/20FRLCD Auto-reader has an internal website that allows knowing in real time and remotely the status, configuration and history of results of the device. You can access the website from any mobile device, tablet or PC connected to the same network as the auto-reader. To access, you must first know the IP address that the router in your network has assigned to the auto-reader. For this, you must press the icon  on the main screen. A screen will be displayed where the batch and serial number of the equipment and its IP address will be informed, besides other data.

## Homepage | Status screen

In any web browser (from a connected device to the same network as the auto-reader) enter the following address: <http://xxx.xxx.xxx.xxx>. The "X" characters represent the IP number reported on the device. Then, the equipment's website homepage will open. On this screen you will be able to see the status of each of the positions as well as the remaining incubation time of those in the middle of a reading process. Also, from this screen, you will have access to the history of results  and settings .

## Results history

You can access this page from the homepage or from any other page, by clicking on the icon . On this page you will find the history of the last 208 results grouped by 10 per page. By clicking on *Previous* or *Next*, you will be able to scroll all the results that are stored in chronological order.

## Settings

You can access this page from the homepage or from any other page, by clicking on the icon . On this page you will find language configuration, time zone and time. These parameters can be configured and synchronized immediately with the auto-reader.

## Troubleshooting chart

Fault	Possible cause	Action
The auto-reader does not start	Power source is not connected.	Check that the power source is connected to the power supply and that the plug is connected to the auto-reader.
The auto-reader gives an error in one position during "Autotest".	A BI is placed in the position at the moment of starting the auto-reader.	Check that every position is empty at the moment of starting the auto-reader.
The auto-reader gives an error in the "Autotest".	Soil particles might be obstructing the light path between the sensors.	Ensure cleanliness of the incubation position. Use of air is recommended. Do not use solid objects that could damage internal components. Restart the auto-reader.
The auto-reader does not run a BI reading. The red light of that position is blinking.	The autotest gives an error in that position, which becomes disabled.	Make sure that the position is empty when starting the auto-reader. Restart the auto-reader.
The auto-reader does not run a BI readout in any position.	The incubation temperature is not stable.	Wait for temperature to be stable.
The auto-reader does not run Chem dye® Pro readings. The red light of that position is blinking.	The autotest gives an error in the PRO position, which becomes disabled.	Make sure that the PRO position is empty when starting the auto-reader. Restart the auto-reader.
When initiating the PRO incubation program the machine makes a sound but it does not initiate.	The incubation temperature is not stable.	Wait for temperature to be stable.
The auto-reader does not allow to change the incubation program.	Ongoing reading.	Wait for the readings to finish.
The printer does not print and the printer indicating light blinks.	The printer cover is not tightly locked.  Printer without paper.	Check that the cover is tightly closed.  Place a new paper roll in the right direction.
The printer releases unprinted paper.	Paper roll is placed in the wrong direction.	Check the paper roll direction.
The Bionova® traceability software does not detect the auto-reader.	The auto-reader is turned off or network cable is disconnected.	Verify that the equipment is turned on and the network cable is correctly connected.
The update of the auto-reader fails.	The auto-reader is connected to the Bionova® traceability software.  The firewall of the PC operating system prevents the correct functioning of the Bioupdate software.  Some antivirus programs block the update of the auto-reader when running the Bioupdate software.  The auto-reader Information screen is not open.	Close the program, restart the auto-reader and try again.  Disable the PC's firewall for a moment. Once the update process is complete, activate it again.  Disable the Antivirus of your PC for a moment. Once the update process is complete, activate it again.  Before starting your update access the auto-reader Information screen from the main screen by pressing the Bionova icon.
	The connection of the equipment to the network to which your PC is connected fails.	Connect the auto-reader to the network through the specified network cable. Avoid connecting your PC to the network through a Wi-Fi connection. A wired connection to the network is recommended. Make sure your internet service is working correctly.

## **Warranty**

Terragene® S.A. guarantees both the quality of the product material components and the quality of its manufacturing process. Should any material or manufacturing faults be detected within the warranty period (1 year from the time of purchase), the only obligation of Terragene® S.A. will be product repairing or substitution.

### **Limitation of liability**

Terragene® S.A. shall not be liable for any loss or damage that may result from the unsuitable use of the equipment, negligence or user's full responsibility.

### **Intellectual property**

Bionova® IC10/20FRLCD Auto-reader may affect the intellectual property of other manufactures in some countries, invading certain aspects of patents related to this product. Terragene® S.A. shall not be liable for any damages or injuries that may result from the marketing of this product in those countries where Terragene® S.A. does not market it directly.

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Terragene S.A.  
Güemes 2879 (S2002JBM)  
Rosario. Santa Fe. Argentina.

# Index

<b>Composition.....</b>	<b>3</b>
References.....	3
Main screen .....	4
Screen References.....	4
Description of use .....	4
<b>Safety information .....</b>	<b>5</b>
Symbols .....	5
Power supply specifications .....	6
Environment operating conditions .....	6
Serial/Batch Number .....	6
Record Table .....	6
<b>Instructions for use.....</b>	<b>7</b>
Start-up .....	7
Thermal paper specifications.....	8
Replacement of the paper roll.....	8
Remaining incubation time   Configured programs .....	9
Positive control.....	9
Disposal.....	9
Interpretation of results .....	9
Hygiene monitoring system.....	9
Results history visualization .....	11
Audible alarm .....	11
Temperature monitoring .....	11
Care and cleaning .....	11
Network cable specifications .....	11
Firmware update .....	11
Calibration.....	12
<b>Website of the equipment .....</b>	<b>12</b>
Homepage   Status screen .....	12
Results history .....	12
Settings .....	12
<b>Troubleshooting chart .....</b>	<b>13</b>
<b>Warranty.....</b>	<b>14</b>
Limitation of liability .....	14
Intellectual property .....	14
Technical assistance .....	14

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