# IC10/20FRLCD

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# IC10/20FRLCD

## Auto-Reader

*For incubation and readout of Rapid, Super Rapid or Ultra Rapid Biological Indicators and Hygiene Monitoring Systems*

## Usage

For incubation and readout of:

Bionova® Rapid Fluorescence Readout Biological Indicators:  
BT220 (Steam): 3 hs. at 60 °C, BT221 (Steam): 1 h. at 60 °C, BT110 (EO): 4 hs. at 37 °C, BT102 (Formaldehyde): 2 hs. at 60°C, BT95 (VH<sub>2</sub>O<sub>2</sub>): 2 hs. at 60°C.

Bionova® Super Rapid Fluorescence Readout Biological Indicators:  
BT222 (Steam): 1 h. at 60 °C, BT223 (Steam): 30 min. at 60 °C, BT96 (VH<sub>2</sub>O<sub>2</sub>): 30 min. at 60°C.

Bionova® Ultra Rapid Fluorescence Readout Biological Indicators:  
BT224 (Steam): 20 min. at 60°C.

Chemdye® Terragene® Protein Pen System: 15 min at 60 °C

## Applicable Regulation

European Union Low Voltage Directive 2014/35/EU.  
Electromagnetic Compatibility Directive 2014/30/EC.  
RoHS Directive 2011/65/EU.

## Authorization

Designed under Quality Management System standards ISO 13485:2003/NS-EN ISO 13485:2012.

## Characteristics

12 metal positions (3.8 cm deep, 0.9 cm diameter) to incubate Rapid/Super Rapid/Ultra Rapid Fluorescence Readout Self-contained biological indicators, 1 position (3.8 cm deep, 1.3 cm diameter) to incubate Protein Pen, 1 ampoule crusher, 1 special hole for external thermometer.

3.5 " touchscreen LCD. Visual information of the remaining time of the incubation process.

Temperature display on the screen and remaining time for each position.

Dual Temperature System: Allows selecting between two different incubation temperatures (37°C y 60°C).

Allows running different incubation times simultaneously, provided they are of equal incubation temperature.

Audible alarm for event indications.

ETHERNET connection for recording results on PC using a Reading and Traceability software.

Remote display of the screen in PC and Smartphones.

It stores the last 208 results.

Thermal Printer of Results. Allows to reprint the last 208 results.

Automatic detection and cancelation of biological indicator fluorescence reading.

Firmware update by the user.

Fluorescence reading wavelength: ~ 460 nm

Dimensions: 18.0 cm high, 26.0 cm diameter.

Weight: 1.92 Kg.

Voltage range: 100 – 240 V

Power: 28 W

Frequency: 50-60 Hz

### **Environmental conditions during manufacture**

Room temperature 15-35 °C, RH 35-85 %

### **Storage conditions**

T = 10-30 °C, RH 30-80 %, keep away from light.

### **Transportation conditions**

Storage conditions should be strictly followed. The product should be transported in closed and reinforced boxes in order to avoid damages. Product transportation does not represent any risk for human health.

### **Shelf-life**

10 years

### **Packing**

1 unit per box.

Packaging information: product code and description, MAC number, and manufacturer information.

### **Labelling**

Same label on product and product's box: product code and description, batch number, serial number, MAC number, and manufacturer information.

### **Possible target markets**

Healthcare, Food, Pharmaceutical and medical devices industries.

### **Other important information (optional)**

The IC10/20FRLCD incubator has an icon on the incubator screen that allows to select the temperature at 37 °C or 60 °C.

In addition, it allows the configuration of each reading position with a different incubation program.

The IC10 / 20FRLCD incubator has an automatic readout system and a printer to record the data.

It also has a history of records to query the last 208 results, which can be viewed from the device screen or remotely accessing the embedded web server.

NOTE: When necessary, limits of each product's parameters should be specified.

### **Description of identified and/or predictable risks which could represent a risk for the product**

Do not pour any liquid inside.

Do not immerse into any liquid.

Only for indoor use.

Disconnect the power cord before cleaning.

Do not use abrasive or corrosive cleaners or disinfectants.

Make sure that the incubator is connected to a properly rated power cord.