IT12

Chemical Integrator

*For Ethylene Oxide sterilization processes.
Two levels.*

Usage

Control of Ethylene Oxide sterilization processes.

Applicable Regulation

Designed under Quality Management Systems standards ISO 13485:2016/NS-EN ISO 13485:2016.
ISO 11140-1:2014

Classification

Class 1, according to risk.

Authorization

ANMAT (Argentinean National Administration of Drugs, Food and Medical Devices) PM 1614-4.

FDA 510(k)

K191021


Characteristics

Type 5 Integrating Indicator (according to ISO 11140-1).
Card printed with indicator ink, laminated on the back side.
Size: 60 mm x 25 mm.

Level 1 – Red: Exposure level, it turns to yellow to indicate the exposure to the process but it does not assure that the conditions indicated in the integrator have been met (UNSAFE condition).

Level 2 – Purple/Brown: Integration level, it turns to green thus indicating that conditions indicated in the integrator have been met (SAFE condition). These conditions are calibrated with the kill time of a 10^6 *Bacillus atrophaeus* spores population, calculated in BIER (Biological Indicator Evaluator Resistometer). Conditions: EO 600 mg/l, RH 60 %, 54 °C.

100 % Toxic Heavy Metal free.

Environmental conditions during manufacture

T= 15-30 °C, RH = 30-80 %.

Storage conditions

T= 10-30 °C, RH = 30-80 %, keep out of direct light.

Transport conditions

Storage conditions should be strictly followed.
Transport in closed and reinforced boxes in order to avoid damages.
The transport of this product does not represent a risk for health.

Shelf life

5 years.

Packing

250 units per aluminum foil bag.

Packing information: product description, instructions for use, regulation, storage conditions, manufacturer information and data on label.

Labelling

On product's packing: product code and description, process for intended use, classification according to regulation, color change, presentation, batch number, manufacture and expiration date, bar code and datamatrix code.

Possible target markets

Healthcare, Food, Pharmaceutical and Medical Industries.

Other important information

Read product's instructions for use thoroughly before use.

Precautions

Do not store the product near sterilizing agents.

Do not expose this product to Steam, Dry Heat or any sterilization processes other than EO.

Do not reuse the sterilizer until the indicator has changed into the color indicated in the product.

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