PCD225-C

R04 / 02.2026

Process Challenge Device (PCD)

Quality certification

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STEAM

LOT

BT225 SCBI

Geobacillus stearothermophilus ATCC® 7953

LOT

Population / Población

CFU / UFC

D-value / Valor D (121 °C) (*)

min.

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Z-value / Valor Z

°C

Parameters determined at time of manufacture according to ISO 11138-1:2017 and ISO 11138-3:2017 standards. The values shown are reproducible only under the same conditions under which they were determined.

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Terragene® and Bionova® are registered trademarks of Terragene S.A. ATCC® is a registered trademark of American Type Culture Collection.

* Determined by the fraction negative method

IT26-C Chemical Integrator

LOT

Stated Values (SV)

Temperature

Time

121 °C

min.

128 °C

min.

135 °C

min.

SV – stated value: Value of a critical process variable at which the indicator is designed to reach its endpoint as defined by the manufacturer.

According to ISO 11140-1:2014 standard

Lic. Adrián J. Rovetto
Technical Director

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CE Mark

CE REP

Authorized Representative in European Community / European Union

MD

Medical Device

LOT

Batch code

LM

Date of manufacture

Use-by date

Manufacturer

ATTENTION: If your chemical indicator has this logo printed on it, it means that your product is compatible with exclusive Trazanto Automatic Quality Control and Traceability System for cleaning, disinfection and sterilization processes. This System interprets the indicator's results in an accurate and automatic manner, and stores all the process data through Bionova® Software. Use Trazanto scanner to interpret the results (see instructions for use in the manual of the device). ATTENTION: Do not write on the indicator before scanning.

Result reference guide

Unprocessed

Accept

Reject

EN

Process Challenge Device

For steam sterilization processes

Indications for use
United States
Terragene® Bionova® Photon® Process Challenge Device with moving-front integrator (PCD225-C) provides a defined challenge resistance against the claimed cycles shown below and demonstrated resistance equivalence to the ANSI/AAMI 16 towel pack. The device provides routine monitoring of the following steam sterilization processes:
• *Gravity-displacement Steam Sterilization Cycles*
132 °C, 25 minutes
132 °C, 15 minutes
132 °C, 10 minutes
135 °C, 10 minutes
• *Dynamic-air-removal Steam Sterilization Cycles*
132 °C, 4 minutes
135 °C, 3 minutes
Terragene® Bionova® Photon® Auto-reader Incubator (BPH) incubates at 60 °C and reads the Terragene® Bionova® Photon® Biological Indicator (BT225) which is a 7-second readout Self-Contained Biological Indicator (SCBI) inoculated with a minimum of 10⁶ viable *Geobacillus stearothermophilus* bacterial spores.
Outside the United States
Terragene® Bionova® Photon® Process Challenge Device with moving-front integrator (PCD225-C) has been designed for quick and easy routine monitoring and qualification testing of dynamic-air-removal and gravity displacement steam sterilization processes at 132-135 °C.

Device description
Terragene® Bionova® Photon® Process Challenge Device with moving-front integrator (PCD225-C) has been designed to simulate a load to be sterilized and to pose a challenge to the sterilization process. It is used to evaluate the effective performance of the process by detecting inadequate air removal and steam penetration. It also allows release of routine loads, especially implants.
Terragene® Bionova® Photon® Process Challenge Device with moving-front integrator (PCD225-C) consists of a disposable pre-assembled test pack as outlined in ANSI/AAMI ST79 which contains a Bionova® BT225 Self-Contained Biological Indicator (SCBI), a PCDBI-C-RC Record Card and an Integron® IT26-C moving front chemical integrator (Type 5 according to ISO 11140-1:2014 standard) that gives instant visible indication that sterilizing conditions have been reached. Each PCD consists of a stack of porous cards holding a SCBI that contains a population of *Geobacillus stearothermophilus* ATCC® 7953 spores soaked on a carrier as well as growth indicator medium contained in a glass ampoule. Each SCBI has a process indicator (Type 1 according to ISO 11140-1:2014 standard) on the label that changes from pink to brown when exposed to steam. The migrating chemical integrator shows ACCEPT result when sterilization conditions were reached while the process indicator (Type 1 according to ISO 11140-1:2014 standard) on the PCD box changes from light blue to dark grey/black when exposed to steam.

Precautions
WARNING: Do not use the PCD for monitoring hydrogen peroxide, ethylene oxide, dry heat, formaldehyde or any sterilization process other than steam. Do not reuse SCBIs.
WARNING: Place one or more PCDs in sterilizing hard-to-reach areas to ensure all areas of the chamber are sterilized. Evaluate all load configurations to ensure all hard-to-reach areas have been identified, and place a PCD in each of those locations.

Instructions for use
1. Place the test pack inside a normally loaded steam autoclave, in those areas which are considered most inaccessible for the sterilizing agent (e.g., the center of the load and areas near the door).
2. Run the sterilization cycle.
3. After the sterilization process has finished, open the sterilizer door, wait for 5 minutes and remove the PCD. **NOTE:** The color of the box may vary from the original after undergoing the sterilization cycle. This does not represent a problem regarding the operation or quality of the product.
4. Check that the process indicator printed on the box has changed color from light blue to dark grey/black. Open the test

pack, wait 5 minutes and remove the SCBI. Allow it to cool down to room temperature. **PRECAUTION:** Wear safety glasses and gloves when removing the SCBI from the sterilized test pack. **WARNING:** Do not crush or handle the SCBI excessively, since this might cause the glass ampoule to burst.
5. Check Integron® IT26-C moving-front chemical integrator for correct exposure. If the dark bar has reached the ACCEPT zone, this confirms that the inside of the pack has been exposed to correct sterilization conditions (refer to Result Reference Guide). Otherwise, check the sterilization process.
6. Check the process indicator printed on the SCBI's label. A color change to brown confirms that the SCBI has been exposed to steam, **IMPORTANT:** This color change does not indicate that the process was sufficient to achieve sterility.
7. Identify the Bionova® BT225 SCBI by writing the sterilizer number (in case of having more than one), load number, and processing date on the label. Fill out the required information on the Record Card.
8. Activate the SCBI by pressing the lid to seal the tube and crushing the ampoule contained in it with an individual ampoule crusher or with the ampoule crusher placed in the rear of the Bionova® Photon® Auto-reader Incubator (BPH). Then shake the tube down vigorously, until the medium reaches the base of the tube and soaks the spore carrier entirely. Incomplete wetting of the spore carrier may lead to an incorrect fluorescence readout. **NOTE:** Holding time between sterilization and activation should not exceed a 7-day period. **IMPORTANT:** Use a non-sterilized SCBI as a positive control at least once per day, when a sterilization cycle is run. The positive control ensures that the correct incubation conditions have been followed, that the medium promotes rapid growth, that spore viability has not been compromised by improper storage temperature, humidity or proximity to chemicals, and that the Bionova® Photon® Auto-reader Incubator (BPH) is functioning properly. Both the positive control indicator and the processed indicator should belong to the same batch.
9. Immediately after the activation, incubate processed and control biological indicators in the Bionova® Photon® Auto-reader Incubator (BPH) for 7 seconds at 60 ± 2 °C to get the final fluorescence result.
WARNING: Bionova® Photon® Auto-reader Incubator (BPH) can detect residual fluorescence signal even when the *G. stearothermophilus* spores are damaged to the point they will not grow. Therefore, a SCBI may be read as fluorescence-positive even when it is visually growth-negative. This characteristic of the fluorescence readout response provides an increased margin of safety and meets the criteria of the FDA reduced incubation time protocol. For this reason, the incidence of false positive results for Bionova® Photon® Biological Indicator (BT225) may be higher compared to a typical SCBI in partial or ineffective sterilization cycles. **WARNING:** In case of a 7-second fluorescence readout positive result after a full sterilization cycle, we recommend reprocessing the load.
Fluorescence detection by the Bionova® Photon® Auto-reader Incubator (BPH) (excitation 340-380 nm / emission 455-465 nm) means a failure in the sterilization process. If no fluorescence is detected at 7-second incubation, the result is negative (i.e., the sterilization process was successful). The positive control must give positive fluorescence readout for results to be valid. **NOTE:** A 48-hour readout is optional to confirm the 7-second result by visual color change. If the sterilization process has not been successful, culture medium will turn yellow during incubation at 60 ± 2 °C. If sterilization was successful, culture medium will remain purple after incubation. It is good practice to incubate a positive control for a visual color change. The positive control must show a color change from purple to yellow for results to be valid.
10. Record the Bionova® Photon® Biological Indicator (BT225) and Integron® IT26-C results and adhere the self-adhesive Record Card or, alternatively, only the area containing the integrator indicator.
11. Discard the test pack and the indicators immediately, as it is indicated below.
NOTE: If any serious incident occurs in relation to the device, it should be reported to Terragene S.A. and the competent authority of the State in which the user is established.

Monitoring frequency
Follow the policies and procedures with the monitoring frequency specified by the professional associations and/or standards corresponding to your country.
Storage
Store in a dark place, at a temperature between 10-30 °C, 30-80 % relative humidity. Do not freeze. Do not store near sterilizing agents or other chemical products.
Shelf life
Bionova® Photon® Process Challenge Device with moving-front integrator (PCD225-C) has an expiration date of 18 months from the date of manufacture, given by the SCBI that carries, when stored at recommended conditions. Do not use it after expiration date. Chemical integrators and process indicators have an expiration date of 18 months when used in/on PCD. Endpoint stability reaction: Chemical indicator endpoint shall remain unchanged for a period of not less than 6 months when stored at previously indicated conditions.

Disposal
Discard biological indicators after use according to your country's healthcare and safety regulations. The positive SCBIs can be autoclaved in a gravity air displacement steam sterilizer at 121 °C for 30 minutes, 132 °C for 15 minutes or 134 °C for 10 minutes; or in a dynamic air removal steam sterilizer at 132 °C for 4 minutes or 135 °C for 3 minutes.
Discard the test pack after use along with paper waste, according to your country's healthcare and safety regulations.

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