
Challenge Kit for Endoscope Reprocessors

Results reference guide

Unprocessed

Failures

- Inadequate chemical agents' concentration.
- Insufficient water pressure.
- Poor water quality, temperature or cycle duration.

Processed and satisfactory

NOTE: Reference colors exhibited in printed prospects as well as those displayed in our website and soft copies of other documents, do not necessarily represent the real color as shown in the actual indicators. Images are reproduced for illustration purposes only, and do not necessarily match the article described in every detail.

**ATTENTION:** If your chemical indicator has this logo printed on it, it means that your product is compatible with exclusive **Trazanto** Automatic Quality Control and Traceability System for cleaning, disinfection and sterilization processes.
This System interprets the indicator's results in an accurate and automatic manner, and stores all the process data through **Bionova® Software**. Use **Trazanto** scanner to interpret the results (see instructions for use in the manual of the device).

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Explanation of Symbols

 Medical Device.
 Batch number.
 Manufacture Date.
 Expiration Date.
 Manufacturer.

 CE mark.
 Authorized representative in the European Community.

Figure 1


Figure 2


Figure 3



Figure 4

Cleaning Challenge Kit

For monitoring cleaning performance of endoscopes in Automated Endoscope Reprocessors

Endoscopes are highly-complex reusable medical devices. They have different channels, up to 2 meters long, with different inner diameters. Since these devices are in contact with human tissues and mucosa, they must undergo cleaning and high-level disinfection processes before being reused, in order to avoid the transmission of pathogens. Automated Endoscope Reprocessors (AERs) offer the cleaning and disinfection method par excellence, since they require minimal user intervention, and, also, have highly reproducible cleaning and disinfection cycles. An optimal result depends on the proper combination of water pressure and quality, an effective concentration of chemical agents, temperature and cycle duration. The monitoring of cleaning performance is essential for guaranteeing compliance of cleaning standards.

Product description

Chemdye® Lumenia Kit consists of a challenge device (Chemdye® Lumenia CCDER) and cleaning indicators (Chemdye® Lumenia CDWE) for monitoring cleaning performance in Automated Endoscope Reprocessors. Chemdye® Lumenia CCDER consists of a device connected to a high strength and durability PTFE (polytetrafluoroethylene) tube that simulates an endoscope's lumen, with length 2 m and an inner diameter of 1 mm or 2 mm, depending on the device model:

- Chemdye® Lumenia CCDER1: a device and a PTFE tube with an inner diameter of 1 mm.
- Chemdye® Lumenia CCDER2: a device and a PTFE tube with an inner diameter of 2 mm.
- Chemdye® Lumenia CCDER122: three devices, a PTFE tube with an inner diameter of 1 mm and two PTFE tubes with an inner diameter of 2 mm.

Each PTFE tube has an input terminal for plugging it to the Automated Endoscope Repressor. Two input terminals are included, one male and one female. Select the most appropriate alternative to the terminal of the washing machine. The device allows to hold two Chemdye® Lumenia CDWE cleaning indicators (one inside and one outside). Chemdye® Lumenia CDWE cleaning indicator consists of a synthetic, self-adhesive and thermostable support of 38 mm x 6 mm, printed with a test soil made up from a mixture of especially combined colored natural components. The test soil was developed and validated in compliance with the requirements described in ISO/TS 15883-5:2005 standard. The cleaning indicator's formulation allows the detection of failures that may affect the cleaning performance of endoscopes (both internally and externally) in Automated Endoscope Reprocessors. The test soil must be completely washed in the internal (A) as well as in the external (B) side of the device (Fig. 1). It is recommended to use both indicators, although optionally, you can verify only the internal cleaning. The indicator's ink is not toxic.

Precautions

- Do not reuse cleaning indicators.
- Do not touch the indicator ink with your hands.

Instructions for use

1. Check that Chemdye® Lumenia CCDER Device is clean, dry and well-functioning.
 2. Check that your hands are clean and dry.
 3. Unscrew the devices' cap.
 4. Place a Chemdye® Lumenia CDWE cleaning indicator in the inner cavity of the device in a way that the test soil is facing upwards (Fig.2). **WARNING:** Before placing the indicator in the device, remove the self-adhesive paper. Make sure that the indicator is firmly adhered to the device by using the provided tweezer (Fig.3).
 5. Screw the cap again.
 6. Place a Chemdye® Lumenia CDWE cleaning indicator in the device's external slot in a way that the ink is visible through its window (Fig.4). **WARNING:** Before placing the indicator in the slot, unlock the device (Fig. 4) and remove the self-adhesive paper. Make sure that the indicator is firmly adhered to the slot by using the provided tweezer. Secure the indicator with the device's lock (Fig. 4).
 7. Connect the Chemdye® Lumenia CCDER Device to the Automated Endoscope Repressor through the output terminals.
 8. Run a normal cleaning cycle and cancel it immediately after the cleaning phase and just before the disinfection phase. If it is not possible to interrupt the cycle, check the results when the complete cycle finishes.
 9. Once the cycle is completed, unplug the connectors and remove the Chemdye® Lumenia CCDER Device from the washing machine.
 10. Carefully remove the indicators from the devices, by using the given tweezer.
 11. Examine the indicators, in order to find remains of test soil by placing them on a white background. If you have the Trazanto traceability software, place the indicator in the scanner and proceed with the results' recording. Record the results and adhere the indicator in the record notebook.
 12. If any of the indicators has test soil remains, take note of the indicators that have failed, and take the applicable measures in case of failure of the cleaning performance test.
- NOTE:** If any serious incident occurs in relation to the device, it should be reported to Terragene S.A. and the competent authority of the State in which the user is established.

Interpretation of results

If the Chemdye® Lumenia CDWE cleaning indicator is completely washed, it means that the cycle conditions are appropriate for guaranteeing a correct endoscope's cleaning (both internally and externally). Therefore, the flexible endoscope can be reprocessed by using the recently validated program. The failures on a cleaning cycle may result from an inadequate concentration of chemical agents, insufficient water pressure, inadequate water quality, wrong temperature or duration of the cycle, among others. These facts decrease the cleaning performance, therefore it is not recommended to reprocess or use the endoscope until a satisfactory result is obtained.

NOTE: Chemdye® Lumenia Kit allows the monitoring of endoscopes' cleaning processes, but it does not make it possible to determine if the appropriate disinfection level has been reached. If your equipment does also perform a high-level disinfection by using peracetic acid as disinfecting agent, you can complete the monitoring through the use of Integron® IT400 and IT401 chemical indicators, as well as through Bionova® BT400 biological indicator.

Additional technical information

You should verify that the connectors and the PTFE tubes are in perfect conditions in each use. In case the connector is damaged and its replacement is justified, please contact a representative and request a spare part. In case the tube is blocked or it has water leaks, you should replace the device with a new one. One or three replacement terminal connectors are included depending on the challenge device model.

WARNING: Terragene S.A. does not guarantee the well-functioning of Chemdye® Lumenia Kit if other manufacturers' cleaning indicators or challenge devices are used. Chemdye® Lumenia results are only valid for a cleaning cycle tested with these particular parameters. Do not draw conclusions for cycles with different parameters. Terragene S.A. does not guarantee the well-functioning of Chemdye® Lumenia Kit if it is used for a cleaning process different from the one it was originally designed for.

Storage

Store the Chemdye® Lumenia Kit in a dry place, protected from light, at a temperature from 10-30 °C and 30-80 % relative humidity. Do not use Chemdye® Lumenia CDWE indicators after their expiration date. **NOTE:** An incorrect storage of the product may affect its stability as well as the obtained results.

Disposal

Discard the kit components after their use, according to your country's healthcare and safety regulations.

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