Infection Control Division | **Accessories**

Bionova® MiniBio Auto-reader Incubator

For incubation and readout of Rapid, Super Rapid or Ultra Rapid Biological Indicators

Usage

United States

Terragene® Bionova® MiniBio Reader Incubators incubate at 60 °C and 37 °C and read the Terragene® Bionova® SCBI for fluorescent results at the times prescribed in the User Manual.

Outside the United States

For incubation and readout of:

Bionova® Rapid Fluorescence Readout Biological Indicators:

-BT220 (steam): 3 h at 60 °C.

-BT110 (ethylene oxide): 4 h at 37 °C.

-BT102 (formaldehyde): 2 h at 60 °C.

Bionova® Super Rapid Fluorescence Readout Biological Indicators:

-BT222 (steam): 1 h at 60 °C.

-BT96 (VH₂O₂): 30 min at 60 °C.

-BT97 (VH₂O₂): 1 h at 60 °C.

Bionova® Ultra Rapid Fluorescence Readout Biological Indicators:

-BT224 (steam): 20 min at 60 °C.

Applicable regulation

Designed under Quality Management System standards
ISO 13485:2016/EN ISO 13485:2016

IEC 61010-1

IEC 61010-2-010

Low Voltage Directive 2014/35/EU

EN IEC 61326-1

EMC Directive 2014/30/EU

RoHS Directive 2011/65/EU

WEEE Directive 2012/19/EU

FDA 510(k)

K191021.

Characteristics

3 positions (38 mm deep, 9 mm diameter) to incubate Rapid, Super Rapid and Ultra Rapid Fluorescence Readout Self-Contained Biological Indicators.

1 ampoule crusher.

Hole for external thermometer.

Built-in top cover.

Dual temperature system: Allows selection between two different incubation temperatures (37 °C and 60 °C).

USB connection for recording results in a PC through Bionova® Cloud Web Application.

Allows incubating 3 SCBIs simultaneously, in different incubation times, at the same temperature.

Thermal printer to record results.

Sound alarm for event indication.

Automatic detection and cancellation of biological indicator fluorescence reading.

Firmware update by the user.

Infection Control Division | **Accessories**

Bionova® MiniBio Auto-reader Incubator

For incubation and readout of Rapid, Super Rapid or Ultra Rapid Biological Indicators

Fluorescence reading wavelength: ~ 460 nm.
Dimensions: 105 mm high, 146 mm wide and 155 mm deep.
Voltage range: 100-240 VAC.
Power: 24 W.
Frequency: 50-60 Hz.
Operating temperature: 10-30 °C.

Storage conditions

T = 5-40 °C, RH = 30-80 %, keep away from direct sunlight.

Transport conditions

It is recommended to respect the storage conditions during transport.

Products should be transported in closed and reinforced boxes in order to avoid damages.

The transport of this product does not represent any risk for human health.

Expected lifetime

5 years.

Packing

1 unit per box.

Packing information: Product code and description, regulatory data, storage conditions and manufacturer information.

Labelling

On product: Product code and description, operating conditions, manufacture date, manufacturer information, batch number, serial number, manufactory firmware and Data Matrix.

On product's box: Product code and description, storage conditions, regulatory information, manufacture date, manufacturer information, batch number, serial number, manufactory firmware and Data Matrix.

Possible target markets

Healthcare, Food, Pharmaceutical and Medical Devices industries.

Other important information

Bionova® MiniBio Auto-reader has a button located on the incubator panel to select temperatures at 37 °C or 60 °C. It has a button for each of the positions of readout to set up a different time for each position.

Infection Control Division | **Accessories**

Bionova® MiniBio Auto-reader Incubator

For incubation and readout of Rapid, Super Rapid or Ultra Rapid Biological Indicators

Precautions

- Do not pour any liquid inside.
- Do not immerse in any liquid.
- For indoor use only.
- Disconnect the power cord before cleaning.
- Do not use abrasive or corrosive cleaners or disinfectants.
- Ensure that the auto-reader is connected to a power cord with the appropriate rating.
- This product contains dry natural rubber.

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