Process Challenge Device (PCD)

Quality certification

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LOT

BT220 SCBI

Geobacillus stearothermophilus ATCC® 7953

LOT

Heat shock population / Población _____ CFU / UFC

D-value / Valor D (121 °C) (*) _____ min.

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Survival time = (log₁₀ labeled population - 2) x labeled D-value

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Kill time = (log₁₀ labeled population + 4) x labeled D-value

Z-value / Valor Z _____ °C

Parameters determined at time of manufacture according to ISO 11138-1:2017, ISO 11138-3:2017 and IRAM 37102:1999 (Parts 1 and 3) standards. The values shown are reproducible only under the same conditions under which they were determined.

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ISO and USP Compliant

Terragene® and Bionova® are registered trademarks of Terragene S.A.

ATCC® is a registered trademark of American Type Culture Collection.

* Determined by the fraction negative method

PCDBI-1-RC Chemical Integrator

LOT

Stated Values (SV)

Temperature

Time

121 °C _____ min.

128 °C _____ min.

135 °C _____ min.

SV – stated value: Value of a critical process variable at which the indicator is designed to reach its endpoint as defined by the manufacturer.

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Explanation of Symbols

STEAM	Product designed for use with Steam sterilization cycles.	MD Medical Device.
CE	CE mark.	LOT Batch number.
EC	Authorized representative in the European Community.	MAN Manufacture Date.
REP		EXP Expiration Date.
		MAF Manufacturer.

EN Process Challenge Device

For a rapid and easy monitoring of Steam sterilization processes

Indications for use

Bionova® PCD220-1 Steam Process Challenge Device has been designed for quick and easy monitoring of dynamic air removal and gravity air displacement Steam sterilization processes at 132/135 °C ≥ 4 minutes and at 121 °C ≥ 30 minutes.

Device description

Bionova® PCD220-1 Process Challenge Device (PCD) Test Pack has been designed to simulate a load to be sterilized and to pose a challenge to the sterilization process. It is used to evaluate the effective performance of the process by detecting inadequate air removal and steam penetration. It also allows release of routine loads, especially implants. Furthermore, it allows to perform the routine monitoring and periodic validation of the sterilizers (after repair, installation, relocation). The device consists of a disposable pre-assembled package as outlined in ANSI/AAMI ST79 and contains a Self-Contained Biological indicator (SCBI) and a Type 5 Chemical integrator indicator printed on The Record Card. The chemical integrator on the Record Card indicates that sterilizing conditions were met. Each pack consists of a stack of porous cards holding a self-contained biological indicator (SCBI) tube that contains a population of *Geobacillus stearothermophilus* ATCC® 7953 spores soaked on a carrier as well as growth indicator medium contained in a glass ampoule. The pack also contains a self-adhesive Record Card with a Type 5 Chemical integrator indicator and, on the outside, a Type 1 process indicator.

Precautions

Do not use the PCD for monitoring ethylene oxide, dry heat, formaldehyde or any sterilization process other than steam. Do not reuse the SCBI.

Instructions for use

- Place the pack inside a steam autoclave with the load to be sterilized. Place it in areas that you consider a priori most inaccessible for the sterilizing agent (steam). Generally, the problematic areas are the center of the load and the bottom shelf near the door and over the drain.
 - Run sterilization cycle.
 - After the sterilization process has finished, open the sterilizer door, wait for 5 minutes and remove the test pack. **NOTE:** The color of the box design may vary from the original color after undergoing the sterilization cycle. This does not represent a problem regarding the operation or quality of the product.
 - Check that the external process indicator color on the outside of the test pack has changed. Open the test pack, wait 5 minutes and remove the SCBI. Allow it to cool down to room temperature. **PRECAUTION:** Wear safety glasses and gloves when removing the biological indicator from the sterilized test pack. **WARNING:** Do not crush or handle the biological indicator excessively, since this might cause the glass ampoule to burst.
 - Check the chemical integrator on the Record Card for correct exposure. Color change to the reference confirms that the inside of the pack has been exposed to sufficient sterilization conditions. Otherwise, check the sterilization process.
 - Check the process indicator on the label of the biological indicator. A color change to brown confirms that the biological indicator has been exposed to steam. **IMPORTANT:** This color change does not indicate that the process was sufficient to achieve sterility.
 - Identify the Bionova® BT220 SCBI by writing the sterilizer number (in case of having more than one sterilizer), load number, and processing date on the label. Fill out the required information on the Record Card.
 - Crush the glass ampoule contained in the SCBI. Then shake the tube down vigorously, with movements similar to those performed to lower the temperature in a mercury thermometer, until the medium reaches the base of the tube and soaks the spores carrier entirely. Incubate at (60 ± 2) °C in an appropriate Bionova® auto-reader incubator. **IMPORTANT:** Use a non-sterilized biological indicator as a positive control whenever a processed biological indicator is incubated. The positive control ensures that appropriate incubation conditions are met. Both the positive control indicator and the processed indicator should belong to the same batch. The indicator used as positive control must yield a positive result.
 - Incubate the processed biological indicator and the indicator used as positive control for 3 hours at (60 ± 2) °C to get the final fluorescence result. A 48-hour readout by visual color change is optional to confirm the 3-hour result. A positive fluorescence result (or growth medium color change after the 48-hour optional incubation) of the processed biological indicator means a sterilization process failure has occurred. If a negative fluorescence result is obtained (or growth medium remains the original color after the optional incubation period) the sterilization process was satisfactory.
 - Wait for the final results. Fill out the SCBI results and adhere the self-adhesive whole Record Card or alternatively, only the area containing the integrator indicator. **WARNING:** Do not use the sterilizer until the biological indicator test results are negative.
 - Discard the test pack and the SCBI immediately.
- NOTE:** If any serious incident occurs in relation to the device, it should be reported to Terragene S.A. and the competent authority of the State in which the user is established.

Monitoring frequency

Follow the policies and procedures with the monitoring frequency specified by the professional associations and/or standards corresponding to your country. As recommended practice, and to provide optimal patient safety, Terragene S.A. recommends that each sterilization load be monitored with the appropriate biological indicator.

Storage

Store in a dark place, at temperatures between 10-30 °C, 30-80 % relative humidity. Do not freeze. Do not store near sterilizing agents or other chemical products.

Shelf Life

PCD220-1 has an expiration date of 2 years from the date of manufacture, given by SCBI that carries, when stored at recommended conditions. Do not use after expiration date. Chemical Integrators and Process Indicators have an expiration date of 2 years when used in/on PCD. Endpoint Stability Reaction: chemical indicator endpoint shall remain unchanged for a period of not less than 6 months when stored at previously indicated conditions.

Disposal

Dispose used biological indicators according to your country's healthcare and safety regulations. The positive biological indicator can be sterilized at 121 °C for 30 minutes, 132 °C for 15 minutes or 134 °C for 10 minutes in a gravity air displacement steam sterilizer, or at 132 °C for 4 minutes or 135 °C for 3 minutes in a dynamic air removal steam sterilizer.

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